

# Bed Raisers – 40mm or 100mm



- Choice of 40mm or 100mm raisers.
- Can be stacked up to 2 high giving 4 different heights:
  - 40mm = 1 x 40mm raiser.
  - 80mm = 2 x 40mm raisers.
  - 100mm = 1 x 100mm raiser.
  - 140mm = 1 x 40mm raiser + 1 x 100mm raiser.
- Suitable for chair and bed legs and castors up to 65mm diameter.
- Combined patient and furniture must be less than 250kg across four raisers.



- Bed Raiser Biscuit – fits inside the well on top of the raisers.
- Can be screwed to the bottom of larger legs so that they are secure on the raisers.
- Bed Raiser Mat – fits inside the well on top of the raisers.
- Reduced chance of legs moving on the raisers.