

ACCORA FloorBed Floorline Beds

FloorBed 1 Plus



FloorBed 2



www.floorbed.com

- Stylish
- Easily Transportable
- Reduced Carer Burden
- Enables greater Risk Management
- Multilevel lockout of handset
- Adjustable mattress guides
- Concealed wheels
- No trailing wires

Selection of mattresses available depending on requirements of user.

	FloorBed 1 Plus	FloorBed 2
Size	King Single – L2280 x W1090mm	Single - L2280 x W940 mm
SWL (Patient Weight)	260kg (225kg)	185kg (150kg)
Back Raise	✓	✓
Knee Break	✓	✓
Trendelenburg & Reverse Tilt	✓	✓
High/Ultra Low	70 – 650 mm height range	70 – 800 mm height range
Available Accessories	<ul style="list-style-type: none"> • Bed Lever • Rotating Bed Lever • Self Help Pole • Side Rails (pair) 	<ul style="list-style-type: none"> • Bed Lever • Rotating Bed Lever • Self Help Pole • Side Rails (pair)
Features	The additional 150mm width gives extra space which may be required when turning.	Reduced head and foot board heights in relation to mattress platform, enabling easy access for the carer.

ACCORA FloorBed Floorline Beds



FloorBed 1 Plus



With Extension Kit
(extends mattress
platform by 20cm)



FloorBed 2

Shown with Side Rails and Bed Lever
(also available for FloorBed 1 Plus)



Rotating Bed Lever



Max. height to facilitate care.



With Extension Kit