



ALIVIO

Masterpiece
Collection

**Lift Recliner
Range**

**STYLE &
COMFORT**



The Alivio Masterpiece Lift Recliner Range are a combination of style, support & comfort.

NEW



Donatello Plus

NEW



Carrington



Picasso



Leonardo



Michaelangelo

Petite



Michaelangelo Petite



Donatello - Lateral



Donatello - Roll



Donatello - Tall

Petite



Donatello Petite - Lateral

Petite



Donatello Petite - Roll



Donatello Canyon

Petite



Donatello Canyon Petite



Kahlo



Introduction

With the Alivio Masterpiece Range of premium Lift Recliners, we've got comfort down to a fine art.

Alivio Lift Recliners sculpt and mould to the contours of your body, providing a range of soft yet supportive positions. Take control of your comfort level with multiple independent support features – including a lift and tilt function for a more secure entry/exit into the recliner.

Like any work of art, the Alivio Range blends modern aesthetics with a range of ultra-durable fabrics – designed to suit any home.

Features

- **Two Stage Lift Mechanism**
has a more vertical lift with a slight tilt compared to the forward tilt you find on generic chairs.
- **Lay-flat Recliner**
- **Footrest Extension**
- **Massage**
– Tri-level massage function
(Donatello Plus and Carrington models only)
- **Hand Controller**
– Easy to use
– Soft grip buttons
– Choose “full recline” or “full stand” with the touch of a button
- **Heated Seat**
– Tri-level heated seat settings
(Donatello Plus model only)
- **Ultra Durable Stylish Fabric**
- **Lumbar Support**
to accommodate sitting posture
- **Power Headrest** to lift or lower the head position
- **Zero Gravity**
(Donatello Plus model only)

General conditions that the power headrest and lumbar adjustment can assist with:



Four* independent motors for a wide range of position combinations and support:

*Kahlo Lift Recliner (KA557) has 3 independent motors. Donatello Plus (KA554-PL) has 5 independent motors.



Power Headrest

The individually adjustable head support allows you to lift or lower the position of the headrest even when the backrest is fully reclined providing a higher level of adjustment for your neck and back.



Power Lumbar

Individually adjustable lumbar support, which can assist in maintaining a proper sitting posture; helping to relieve tension on the spine to prevent pain or discomfort due to poor posture, arthritis or existing back problems.




Power Backrest and Power Footrest

Independently adjust the backrest and the footrest for the position that best fits the clinical needs of the client.



Ultra Durable Fabrics



Easy to Use Controller*

*4 motor controller shown



Range to suit different body types and conditions



Full Lift



Made for comfort



Footrest Extension



Expert Advice



Wireless Phone Charger*

* Donatello Plus and Carrington models only



Soft-Touch Controller

Customer & Therapist Feedback

Customer: "I came over to your shop yesterday to try out the latest recliner to hit Australia and wow what a recliner! There is nothing in your shop that comes close and I've tried them all.

Having a fused neck, arthritis and just had two lower back surgeries I am needing support while sitting. I could never get a recliner to support my lower back and neck until now because of the electric adjustable lumbar support and electric adjustable head rest.

Gone are the days of trying to use pillows to get the support which could never be achieved. Even the way the recliner lifts you has improved as you don't feel like you are

going to land on your face."

Therapist: "For clients with postural complaints or conditions - the ability to customise head, lumbar, back and feet position; can make a world of difference to their comfort levels. The ability to make frequent posture adjustments over a period of time is also very reassuring.

The Alivio Lift Recliner Range provides an excellent choice of sizes, features, and individual preferences - which is why I recommend Alivio to my clients time and time again. "

Donatello Plus Lift Recliner KA554-PL

NEW




Full Lift



Massage



Heated seat



Drinks Cup holder



Heat and Massage settings



Zero Gravity

Fabric Name: Canyon Steel

Specifications

Overall Length 1830mm
Overall Width 1010mm
Overall Depth 1060mm
Overall Height 1090mm
Seat Width 550mm

Seat Depth 530mm
Backrest height 690mm
Floor To Seat Height 495mm
Product Weight 79kg
Max User Weight 158kg

5 independent motors

- Neck
- Lumbar
- Leg
- Back
- Zero gravity

Carrington Lift Recliner

KA558

NEW



Full Lift



Massage



Wireless Charger & Cup Holder



Neck Support

Lumbar Support

Leg Rest Support

Back Support



Soft Touch controller

Fabric Name: Canyon Silt

Specifications

Overall Length	1830mm
Overall Width	1010mm
Overall Depth	1060mm
Overall Height	1090mm
Seat Width	550mm

Seat Depth	530mm
Floor To Seat Height	508mm
Product Weight	85kg
Max User Weight	158kg

4 independent motors

- Neck
- Lumbar
- Leg
- Back

Picasso Lift Recliner

KA551



Full Lift



Upright



Partial Recline



Full Recline

Fabric Name: Steel

Specifications

Overall Length	1830mm
Overall Width	813mm
Overall Depth	1003mm
Overall Height	1080mm
Seat Width	521mm

Seat Depth	508mm
Max Seat Height	495mm
Floor To Seat Height	495mm
Product Weight	74kg
Max User Weight	158kg

4 independent motors

- Neck
- Lumbar
- Leg
- Back

Leonardo Lift Recliner

KA552



Full Lift



Upright



Partial Recline



Full Recline

Fabric Name: Santabarbara Fog

Specifications

Overall Length	1830mm
Overall Width	902mm
Overall Depth	1067mm
Overall Height	1092mm
Seat Width	571mm

Seat Depth	521mm
Max Seat Height	483mm
Floor To Seat Height	483mm
Product Weight	76kg
Max User Weight	158kg

4 independent motors

- Neck
- Lumbar
- Leg
- Back

Michelangelo Lift Recliner

KA553



Full Lift



Upright



Partial Recline



Full Recline

Fabric Name: Badlands Saddle

Specifications

Overall Length 1830mm
Overall Width 890mm
Overall Depth 1080mm
Overall Height 1067mm
Seat Width 559mm

Seat Depth 546mm
Max Seat Height 470mm
Floor To Seat Height 470mm
Product Weight 70kg
Max User Weight 158kg

4 independent motors

- Neck
- Lumbar
- Leg
- Back

Michelangelo Petite Lift Recliner

KA556



Full Lift



Upright



Partial Recline



Full Recline

Fabric Name: Badlands Saddle

Specifications

Overall Length	1830mm
Overall Width	851mm
Overall Depth	978mm
Overall Height	1041mm
Seat Width	457mm

Seat Depth	457mm
Max Seat Height	457mm
Floor To Seat Height	457mm
Product Weight	70kg
Max User Weight	136kg

4 independent motors

- Neck
- Lumbar
- Leg
- Back

Donatello Lift Recliner

KA554



Full Lift



Upright



Partial Recline

The **Donatello Lift Recliner** has a range of configurable back rest options to suit you:



Standard Support
Backrest: KA554-SB



High Support
Backrest: KA554-TB



Lateral Support
Backrest: KA554-LS



Full Recline

Fabric Name: Merino Steel

Specifications

Overall Length	1830mm
Overall Width	890mm
Overall Depth	1016mm
Overall Height	1092mm
Seat Width	520mm

Seat Depth	520mm
Max Seat Height	508mm
Floor To Seat Height	508mm
Product Weight	74kg
Max User Weight	158kg

4 independent motors

- Neck
- Lumbar
- Leg
- Back

Donatello **Petite** Lift Recliner

KA555



Full Lift

Upright

Partial Recline

The **Donatello Petite Lift Recliner** has a range of configurable back rest options to suit you:



Standard Support
Backrest: KA555-SB

Lateral Support
Backrest: KA555-LS



Full Recline

Fabric Name: Merino Steel

Specifications

Overall Length	1830mm
Overall Width	851mm
Overall Depth	978mm
Overall Height	1041mm
Seat Width	457mm

Seat Depth	457mm
Max Seat Height	457mm
Floor To Seat Height	457mm
Product Weight	70kg
Max User Weight	158kg

4 independent motors

- Neck
- Lumbar
- Leg
- Back

Donatello Canyon Lift Recliner

KA554LSC



Full Lift



Upright



Partial Recline



Full Recline

Fabric Name: Canyon Steel

Specifications

Overall Length 1830mm
Overall Width 890mm
Overall Depth 1016mm
Overall Height 1092mm
Seat Width 520mm

Seat Depth 520mm
Max Seat Height 508mm
Floor To Seat Height 508mm
Product Weight 74kg
Max User Weight 158kg

4 independent motors

- Neck
- Lumbar
- Leg
- Back

Donatello **Petite** Canyon Lift Recliner

KA555LSC



Full Lift



Upright



Partial Recline



Full Recline

Fabric Name: Canyon Steel

Specifications

Overall Length	1830mm
Overall Width	851mm
Overall Depth	978mm
Overall Height	1041mm
Seat Width	457mm

Seat Depth	457mm
Max Seat Height	457mm
Floor To Seat Height	457mm
Product Weight	70kg
Max User Weight	158kg

4 independent motors

- Neck
- Lumbar
- Leg
- Back

Kahlo Lift Recliner

KA557



Full Lift



Upright



Partial Recline



Full Recline

Fabric Name: Stonewash Fudge

Specifications

Overall Length	1830mm
Overall Width	890mm
Overall Depth	1016mm
Overall Height	1092mm
Seat Width	520mm

Seat Depth	520mm
Max Seat Height	508mm
Floor To Seat Height	508mm
Product Weight	74kg
Max User Weight	158kg

3 independent motors

- Neck
- Leg
- Back

Improving your life everyday™

At Your Service

Sales, administration
and general enquiries

National 1300 783 783

Your Customer Service Team

customerservice@kcare.com.au

[kcare.com.au](https://www.kcare.com.au)

All information is considered to be true and correct at the date of publication. Changes in circumstances after the time of publication may impact on the accuracy of the information.

Published: June 2023