


# Alternating Air Mattresses

## Softform Premier Active 2:

- Deep contoured centre foam insert.
- Castellated top surface for pressure care and firm edges for stability when standing up / sitting on bed .
- Alternating air cells under foam.
- Suitable for High Risk / Very High Risk when used with an adjustable bed.
- Sizes: Single.
- SWL 247Kg.



## Oska Alto:

- Full mattress replacement system with 20 air cells.
- Static cells in head area to minimise head movement.
- 2-Part cells with alternating top and static bottom section that will never bottom out.
- Pump can be set to static or alternating mode.
- Suitable for Very High Risk when used with an adjustable bed.
- Sizes: Single.
- SWL 175Kg.

