



aspire

LIFT RECLINE CHAIRS

Contemporary luxury & comfort



RAPHAEL QUATTRO LIFT RECLINE CHAIR

Drawing inspiration from the famous Italian artist, the Raphael is the epitome of style & luxury. Infinite positioning can be achieved thanks to its 4-motor engineering; while luxurious microsuede and vinyl upholstery options bring style & class to any environment. A true centre-piece in the modern home.

- Four motor lift recline chair for infinite positioning
- Microsuede and vinyl upholstery options
- Backlit hand control with integrated USB charging point
- Independent adjustable headrest, back, leg & lumbar support functions
- Integrated lift & recline functions for easy sit-to-stand & maximum comfort
- Intuitive headrest support, ideal when watching TV when reclined
- Head & armrest covers for improved longevity
- Available in four sizes: Small, Medium, Large and MAXI
> See page 18 for measurements



* Only available in medium and large sizes.

DA VINCI QUATTRO LIFT RECLINE CHAIR

Da Vinci always captured the crowds & remains the centre-piece in any room. The ultimate in contemporary living, the Da Vinci Quattro is a visually aesthetic offering within the Aspire GOLD lift recline chair range. Finished with soft microsuede upholstery & featuring an intuitive four motor mechanism, style your new lift recline chair with a choice of four colour options to create a relaxation hub in the home.

- Four motor lift recline chair for infinite positioning
- Microsuede modern look
- Backlit hand control with integrated USB charging point
- Independent adjustable headrest, back, leg & lumbar support functions

- Integrated lift & recline functions for easy sit-to-stand & maximum comfort
- Intuitive headrest support, ideal for watching TV when reclined



Pictured in Gunmetal Microsuede



Upright Position



Lift to Stand



Power Lumbar



Back Recline



Leg Recline



Full Recline



Power Headrest

SIGNATURE SPACE SAVER LIFT RECLINE CHAIR

Tailored with premium Warwick Global Suede & Vinyl upholstery & high grade durable foams to provide years of postural support, comfort & safety. The Signature Lift Recline Chair with space saver mechanism is ideal for smaller rooms or confined spaces, conserving the best use of your living space. Available in multiple sizes & three colour options, this lift recline chair is sure to be the *signature* addition to any living space.

- Single motor lift recline chair
- Space saving design mechanism, ideal for small & confined spaces
- Depth adjustable seat
- Backlit hand control with integrated USB charging point
- Available in three sizes, & Warwick Global Suede or Vinyl
**See page 18 for measurements*

COLOUR RANGE



Mink (Microsuede)



Lagoon (Microsuede)



Heritage (Vinyl)



Available in Heritage Vinyl for a
Leather-look feel



Upright Position



Lift to Stand



Recline

Pictured with optional Headrest Cover,
Optional Armrest Covers also available

Relax in a chair that fits

COLOUR RANGE



Mink (Microsuede)



Lagoon (Microsuede)



Heritage (Vinyl)



Pictured with optional Headrest Cover, Optional Armrest Covers also available

SIGNATURE DUAL ACTION LIFT RECLINE CHAIR

The Signature Dual Action Lift Recline Chair offers more positioning, for increased comfort. With independent control of back & legrest functionality, the user has multiple positioning options. Tailored with premium Warwick Global suede & vinyl upholstery & high grade durable foams to provide years of postural support, comfort & safety. Available in multiple sizes & three colour options, this lift recline chair is sure to be the *signature* addition to any living space.

- Dual motor lift recline chair
- Independent control of back & legrest functions
- Backlit hand control with integrated USB charging point
- Depth adjustable seat
- Adjustable backrest padding to support individual posture needs
- Available in three sizes, & Warwick Global Suede or Vinyl
**See page 18 for measurements*



Upright Position



Lift to Stand



Back Recline



Leg Recline



Full Recline

PORTO ZERO GRAVITY PETITE LIFT RECLINE CHAIR

Designed to utilise limited space when fully reclined, without compromising on the overall seat size.

Tilt in Space

NEW Zero Gravity function redistributes sitting pressure, assists with oedema and helps to relax lower back muscles.

- Designed for more petite users
- Head & armrest covers for improved longevity



ZERO GRAVITY



Upright Position



Lift to Stand



Recline



COLOUR RANGE



Fawn (Microsuede)

VERONA SPACE SAVER LIFT RECLINE CHAIR

Bronze RANGE

Comfort you can rely on


Verona signifies the classic Italian lifestyle, encapsulating a classic approach with a modern look & feel. Finished with soft quilted upholstery, this chair offers comfort & style within the home.



COLOUR RANGE



Mushroom (Quilted Upholstery)



MEMORY FOAM

- Single motor lift recliner chair
- Space saving design mechanism, ideal for smaller & confined spaces
- Head & armrest covers for improved longevity
- Memory Foam Seating Surface: Improves comfort and helps to facilitate pressure redistribution



Upright Position



Lift to Stand




Recline



MAUI DUAL ACTION LIFT RECLINE CHAIR

Named after the second largest of the Hawaiian Islands, the Maui Dual Action Recline Chair represents superior beauty & a relaxing lifestyle. A luxurious Lift Recline Chair, the Maui Dual Action delivers multiple positioning options, with additional backrest padding for a plush sitting & reclined experience. Finished with soft quilted upholstery, this chair will be the relaxing staple piece in any home.

- Dual motor lift recline chair
- Independent control of back & legrest functions
- Head & armrest covers for improved longevity
- Memory Foam Seating Surface: Improves comfort and helps to facilitate pressure redistribution



MEMORY FOAM



Upright Position



Lift to Stand



Back Recline



Leg Recline



Full Recline

Bronze RANGE

Comfort you can rely on

COLOUR RANGE



Mocha (Quilted Upholstery)



MONTANA SINGLE ACTION MAXI LIFT RECLINE CHAIR

Bronze RANGE

Comfort you can rely on

Commonly referred to as "Big Sky Country", the sprawling mountain state of Montana epitomises vast space & freedom. The Montana Lift Recline Chair provides a large seating surface, offering improved comfort for users who need that little bit of extra room.

- Single motor lift recline chair
- Designed for larger users
- Head & armrest covers for improved longevity
- Superior wipe clean microsuede upholstery

COLOUR RANGE



Bronze (Microsuede)



Upright Position



Lift to Stand



Full Recline



RANGE OVERVIEW

	GOLD				SILVER			BRONZE					
	RAPHAEL QUATTRO LIFT RECLINE CHAIR		DA VINCI QUATTRO LIFT RECLINE CHAIR		SIGNATURE SPACE SAVER LIFT RECLINE CHAIR	SIGNATURE DUAL ACTION LIFT RECLINE CHAIR		PORTO ZERO GRAVITY PETITE LIFT RECLINE CHAIR	VERONA SPACE SAVER LIFT RECLINE CHAIR	MAUI DUAL ACTION LIFT RECLINE CHAIR	MONTANA SINGLE ACTION MAXI LIFT RECLINE CHAIR		
Description	4 Motor - Quattro				4 Motor - Quattro		1 Motor - Space Saver	2 Motor - Dual Action		Zero Gravity	1 Motor - Space Saver	2 Motor - Dual Action	1 Motor - Heavy Duty
Functions													
Size	Small	Medium	Large	Maxi	Medium		Small	Medium	Large	Petite	Medium	Medium	Maxi
Upholstery	Microsuede	Microsuede & Vinyl		Microsuede	Microsuede		Warwick Global Suede & Vinyl			Soft Suede	Quilted Upholstery	Quilted Upholstery	Microsuede Leather-look
Seat Height	470 mm	470 mm	500 mm	510 mm	480 mm		480 mm	480 mm	500 mm	490 mm	500 mm	500 mm	510 mm
Seat Depth	480 mm	500 mm	530 mm	540 mm	500 mm		460 - 485 mm	490 - 515 mm	510 - 535 mm	500 mm	500 mm	500 mm	540 mm
Seat Width	480 mm	520 mm	550 mm	650 mm	530 mm		480 mm	510 mm	540 mm	490 mm	540 mm	540 mm	650 mm
Overall Height	1060 mm	1090 mm	1080 mm	1050 mm	1050 mm		980 mm	1050 mm	1110 mm	1120 mm	1100 mm	1100 mm	1050 mm
Overall Width	740 mm	780 mm	810 mm	930 mm	870 mm		740 mm	780 mm	780 mm	770 mm	780 mm	780 mm	930 mm
Backrest Height	630 mm	680 mm	700 mm	650 mm	630 mm		580 mm	620 mm	650 mm	690 mm	630 mm	630 mm	650 mm
SWL	130 kg			220 kg	159 kg		145 kg	145 kg	155 kg	150 kg	130 kg	130 kg	220 kg
Warranty	2 Years				2-5 Years *		2 Years			2 Years	2 Years	2 Years	2 Years
AUST R ARTG No.	306075				216484		230186			306075	306075	306075	306075

COLOUR RANGE

- Bronze**
S - CHP227290
M - CHP227300
L - CHP227310
MAXI - CHP227320
- Marina Grey**
S - CHP227291
M - CHP227301
L - CHP227311
MAXI - CHP227321
- Nordic Stone**
M - CHP227302
L - CHP227312

- Latte**
CHP227500
- Fog**
CHP227520
- Dove**
CHP227530
- Gunmetal**
CHP227510

SIGNATURE SPACE SAVER

- Mink**
S - CHP197700
M - CHP197710
L - CHP197720
- Lagoon**
S - CHP197705
M - CHP197715
L - CHP197725
- Heritage - Vinyl**
M - CHP197716
L - CHP197730

SIGNATURE DUAL ACTION

- Mink**
S - CHP197735
M - CHP197745
L - CHP197755
- Lagoon**
S - CHP197740
M - CHP197750
L - CHP197760
- Heritage - Vinyl**
M - CHP197751
L - CHP197765

- Fawn**
CHP227200

- Mushroom**
CHP227100


- Mocha**
CHP227110

- Bronze**
CHP227400

ASPIRE ACCESSORIES

Chair Protector

Universal Size
CHP227600



Signature Chairs - Head & Armrest Cover Set

Mink
CHP197770

Lagoon
CHP197771

Heritage
CHP197772



* 2 Years Handset/Upholstery & Electronics, 5 Years Frame



Contemporary luxury & comfort



ASBR26.44 021221

ASPIRE DEALER

Visit www.aidacare.com.au
or Call 1300 133 120
to Find Your Closest Showroom

