

BC₂

The BC2 day chair is extremely durable and comfortable featuring an ergonomic additional segmented backrest and seat pad with features to provide added comfort and functionality for its users.

Perfectly suited across the aged care, health care and hospital sectors.



Features

Fully welded powdercoated tubular steel frame with height adjustable legs with x 6 height levels, standard finish in light beige or other colours available on request (extended lead times may apply).

Ergonomic additional segmented backrest and seat pad for added comfort and ease of cleaning.

Rear transport wheels and removable meal tray available as added extras.

Upholstered to specification.

Weight capacity of up to 160 Kg.

10 year commercial frame warranty.

Dimensions

Height 990 - 1180 mm **Seat Height Arm Height** 550 - 690 mm 390 - 550 mm Width **Seat Width Arm Width** 100 mm 730 mm 510 mm Depth **Seat Depth** 740 mm 490 mm



BC2 Premium

The BC2 Premium day chair is extremely durable and comfortable featuring an ergonomic additional segmented backrest and seat pad with features to provide added comfort and functionality for its users.

Perfectly suited across the aged care, health care and hospital sectors.



Fully welded powdercoated tubular steel frame with height adjustable legs with x 6 height levels, standard finish in light beige or other colours available on request (extended lead times may apply).

Ergonomic additional segmented backrest and seat pad for added comfort and ease of cleaning.

Wing back support and enclosed arms.

Rear transport wheels and removable meal tray available as added extras.

Upholstered to specification.

Weight capacity of up to 160 Kg.

10 year commercial frame warranty.



Dimensions

Height 990 - 1180 mm **Seat Height** 390 - 550 mm **Arm Height** 550 - 690 mm Width **Seat Width Arm Width** 730 mm 510 mm 100 mm Depth **Seat Depth** 740 mm 490 mm



BC₂ Mobile

The BC2 mobile day chair is extremely durable and comfortable featuring an ergonomic additional segmented backrest and seat pad with features to provide added comfort and functionality for its users.

Perfectly suited across the aged care, health care and hospital sectors.



Fully welded powdercoated tubular steel frame with height adjustable legs with x 6 height levels, standard finish in light beige or other colours available on request (extended lead times may apply).

Ergonomic additional segmented backrest and seat pad for added comfort and ease of cleaning.

100 mm rear locking mobility castors and back handle for ease of movability.

Features slide out mobility footrest and optional meals tray attachment with built in storage for when not in use.

Upholstered to specification.

Weight capacity of up to 150 Kg.

10 year commercial frame warranty.







Dimensions

Height 990 - 1140 mm **Seat Height** 380 - 530 mm **Arm Height** 570 - 710 mm **Seat Width Arm Width** Width 730 mm 510 mm 100 mm Depth 740 mm **Seat Depth** 490 mm



BC₂ King

The BC2 King size day chair is extremely durable and comfortable featuring an ergonomic additional segmented backrest and seat pad with features to provide added comfort and functionality for its users.

Perfectly suited across the aged care, health care and hospital sectors.



Features

Fully welded powdercoated tubular steel frame with height adjustable legs with x 6 height levels, standard finish in light beige, other colours available on request (extended lead times may apply).

Ergonomic additional segmented backrest and seat pad for added comfort and ease of cleaning.

Rear transport wheels and removable meal tray available as added extras.

Upholstered to specification.

Weight capacity of up to 200 Kg.

10 year commercial frame warranty.

Dimensions

Height **Seat Height** 360 - 510 mm **Arm Height** 550 - 690 mm 990 - 1140 mm **Arm Width** Width **Seat Width** 820 mm 610 mm 100 mm Depth **Seat Depth** 740 mm 490 mm



BC2 Super King

The BC2 Super King size day chair is extremely durable and comfortable featuring an ergonomic additional segmented backrest and seat pad with features to provide added comfort and functionality for its users.

Perfectly suited across the aged care, health care and hospital sectors.



Features

Fully welded powdercoated tubular steel frame with height adjustable legs with x 6 height levels, standard finish in light beige, other colours available on request (extended lead times may apply).

Ergonomic additional segmented backrest and seat pad for added comfort and ease of cleaning.

Rear transport wheels and removable meal tray available as added extras.

Upholstered to specification.

Weight capacity of up to 300 Kg.

10 year commercial frame warranty.

Dimensions

Height 990 - 1140 mm **Seat Height** 360 - 510 mm **Arm Height** 550 - 690 mm **Arm Width** Width 980 mm **Seat Width** 100 mm 740 mm **Seat Depth** Depth 740 mm 490 mm



BC Footstool

The BC footstool is a compact and versatile footstool perfectly suited across the aged care, health care and hospital sectors.

Featuring a fully welded powdercoated tubular steel frame, soft cushions top for superior comfort and height adjustable legs.



Features

Featuring a fully welded powdercoated tubular steel frame and height adjustable legs with x 5 height levels, standard finish in light beige, other colours available on request (extended lead times may apply).

Soft cushioned top for superior comfort.

Upholstered to specification.

10 year commercial frame warranty.

Dimensions

Height 340 - 460 mm

Width 420 mm Depth 320 mm