

Atama BC2 Day Chair

BC2 Day Chair

- The BC2 day chair features a 3 segment back cushion which can be adjusted to suit individual requirements.
- The frame is fully welded, and the leg height is adjustable.
- Standard fabrics generally available within 7 days; custom fabrics approx. 12 weeks.
- Warranty: 5 years on metal frame.



Standard Dimensions		King Size Dimensions	
Seat Width	525mm	Seat Width	605mm
Seat Height (minimum)	430mm	Seat Height (minimum)	410mm
Seat Height (maximum)	580mm	Seat Height (maximum)	560mm
SWL	160kg	SWL	200kg