

Decpac Portable Ramps

Model	Max Height Capacity	Length	Width	Weight
Personal 070	11.5cm	70cm	74.5cm	3.9kg
Multipurpose 087	14.5cm	87cm	83.5cm	5.8kg
Multipurpose 100	16.5cm	100cm	83cm	7.5kg
Senior 090	15.0cm	90cm	75cm	6.4kg
Senior 120	19.5cm	120cm	75cm	8.0kg
Senior 135	22.0cm	135cm	75cm	8.1kg
Senior 150	24.5cm	150cm	75cm	10.0kg
Senior 165	13.5cm	165cm	75cm	11.4kg
Senior 185	15.5cm	185cm	75cm	11.9kg
Senior 200	16.5cm	200cm	75cm	12.3kg
Senior 250	20.5cm	250cm	75cm	17.3kg
Senior 300	25.0cm	300cm	75cm	22.1kg
Senior 350	29.0cm	350cm	75cm	24.7kg
EBL 090	15.0cm	90cm	75cm	7.4kg
EBL 120	19.5cm	120cm	75cm	9.7kg
EBL 135	22.0cm	135cm	82cm	11.1kg
EBL 150	24.5cm	150cm	82cm	12.0kg
EBL 165	13.5cm	165cm	82cm	13.3kg
EBL 200	16.5cm	200cm	82cm	16.3kg
EBL 250	20.5cm	250cm	82cm	22.2kg
EBL 300	25.0cm	300cm	82cm	27.3kg
EBL 350	29.0cm	350cm	82cm	28.6kg

Personal – Portable, lightweight ramp. Ideal for climbing kerbs and similar obstacles. Includes carry bag as standard.

Senior – Flexible range of ramps to cope with a range of obstacles. Available in 10 sizes.

Multi-purpose – Longer and wider ramps to suit a wider range of obstacles.

Edge Barrier Limiter (EBL) – Includes an edge barrier on the ramp for safety. Available in 9 sizes.