

# Freedom Folding Electric Wheelchair



- Lightweight folding electric wheelchair - choice of six models.
- Easily folded for transporting. Some models split in two for lifting.
- Suspension and cushioning for a comfortable ride.
- Joystick can be positioned on right, left or at rear (for attendant).
- Single battery pack included. Can be expanded to either two or three battery packs for additional range.

Model	A07 Lite	A06 Classic	A08 Premium	A08L Premium Sport	DE08 Premium Lite	DE08L Premium Lite Sport
SWL	110 kg	115 kg	160 kg	160 kg	160 kg	160 kg
Weight (No Batteries)	20.0 kg	21.5 kg	23.5 kg	24.0 kg	26.8 kg	27.3 kg
Weight - components	n/a	n/a	n/a	n/a	14.3 kg Frame 12.5 kg Drive	14.3 kg Frame 13.0 Drive
Drive Wheel	10in PU	10in PU	10in PU	12in Pneumatic	10in PU	12in Pneumatic
Size unfolded (L x W x H)	95cm x 56cm x 91cm	96cm x 59cm x 93cm	102.5cm x 59cm x 96cm	102cm x 62cm x 90cm	100cm x 60.3cm x 89cm	100cm x 60.3cm x 89cm
Size – folded (L x W x H)	56cm x 65cm x 30cm	59cm x 69cm x 34cm	60cm x 67cm x 32cm	64cm x 74cm x 33cm	60cm x 67cm x 32cm	60cm x 67cm x 32cm
Seat Size (W x D)	41cm x 43cm	44cm x 43cm	45cm x 44cm	49cm x 44cm	49cm x 44cm	49cm x 44cm
Maximum Incline	8 degrees	8 degrees	12 degrees	12 degrees	12 degrees	12 degrees