

Glide Series 2 Leisure Plus Manual Wheelchair



Standard Seating Choices Available

Width: 37cm, 40cm, 43cm, 46cm, 49cm, 52cm, 55cm.

Depth: 37cm, 40cm, 43cm, 46cm, 49cm, 52cm.

- Lightweight manual wheelchair.
- Choice of seat depth and width available - 37cm to 55cm.
- Leg rests - adjustable, detachable, swing away with flip up plastic footplates (as shown, several options available).
- Armrest - full length with moulded arm pads.
- Wheels 24" Diameter rear wheels with 1 3/8inch tyres.
- Castor - 7" MCP (solid).
- Padded backrest Sail Cloth (note: cushion shown is extra).
- 51cm Seat Height standard, can be customised.
- Safe working load 150Kg.
- Choice of colour finishes.

Options and Accessories

- Quick release wheels (pair)
- Cushion, standard.
- Cushion, dartex cover.
- Desk arm (in lieu of full length). Each.
- Adjustable height arm rest. Each.
- Adjustable height desk arm rest. Each.
- Arm pad – Large. Each.
- Trough arm (fixed bracket). Each.
- Trough arm (swivel bracket). Each.
- Seat belt – car type buckle.
- Tray – flip up, fold away.
- Tray – slide in.
- Anti tip bar and rollers (pair).
- Rear tilt bar.
- Oxygen bottle carrier.
- IV pole
- One-arm drive kit (left or right).
- Elevating leg rests with calf pad. Each.
- Vertical leg rests. Each.
- Stump support, 23, 30 or 40cm. Each.
- Heel loops on footplate. Each.
- Aluminium footplates.
- Footplate – one piece, custom.
- Footplate – one piece, flip up.
- Footplate – angle adjustable. Each.
- Footplate no-gap adapter, aluminium.
- External legrest release lever. Each.
- 24 Inch Solid tyres (flat free). Each.
- 24 Inch Aero wheel (6 spoke, black).
- 24 Inch Spider wheel. Each.
- 24 Inch Sports, high pressure. Each.
- 24 Inch BMX. Fixed axle. Each.
- 24 Inch BMX. Quick release. Each.
- 12 Inch Transit wheel kit (pair).