

ICARE IC333

The IC333 Bed is unmatched for function and durability, taking the hospital look out of your room. This four function bed, features dual lift motors totaling 12,000 newtons lifting power, unlimited under-bed clearance and eight heavy duty castors.



KEY FUNCTIONS




Head Lift



Knee Break




Height Adjustable
22-66cm




Trendelenburg and
Reverse Trendelenburg


BED CODES & SIZES



LONG SINGLE
Code: IC333LS
Mattress Surface Size: 2030 x 900
Base Size: 2150 x 960



KING SINGLE
Code: IC333KS
Mattress Surface Size: 2030 x 1070
Base Size: 2150 x 1120




LONG DOUBLE
Code: IC333D
Mattress Surface Size: 2030 x 1350
Base Size: 2150 x 1400

- One hand control




ONE PIECE QUEEN
Code: IC333Q
Mattress Surface Size: 2030 x 1520
Base Size: 2150 x 1580

- One hand control



**SPLIT QUEEN (HEAD & KNEE FUNCTIONS)
(ONE PIECE HI-LO & TRENDELENBERG)**
Code: IC333SQ
Mattress Surface Size: 2030 x 750 (PER SIDE)
Base Size: 2150 x 1580

- Individual hand control per side
- Fits a one piece queen headboard




SPLIT KING
Code: IC333K
Mattress Surface Size: 2030 x 900
Base Size: 2150 x 1920 (as partner)


COLOUR CHOICES



ONYX



STONE



SMOKE



SILVER



VINYL EBONY



VINYL DOVE

NOTE: Always remember that TWO Long Single beds together make a King size bed.

KEY FEATURES



200kg SWL on Long Single & King Single



Strong Breathable Mattress Platform



Available with Headboards & Footboards



External Rechargeable Battery Back-up



Lockable Heavy Duty Castors



300kg SWL on Double and Queen



End and Corner Mattress Retainer Bars




10 Year Warranty



Back Lit Wipeable Hand Piece



No gap between partner mattresses



Unlimited Under Bed Clearance



Ultra Low Design Reducing Risk of Falls



Optional 3 sided surround for wheelchair transfers

TECHNICAL SPECIFICATIONS

- Safe working load - 200kg on long single and king single sizes and 300kg for double and queen sizes. The safe working load includes any model of Icare mattresses.
- Hand Piece - Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed. Wipeable.
- Fabric - Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties.
- Mattress Surface - Strong metal mesh under mattress allowing increased airflow
- Actuators - 1x 6000 newton - 24 watt head adjust actuator. 1x 6000 newton - 24 watt knee adjust actuator. 2x 6000 newton - 24 watt height adjust actuators (2x 8000 newton on double and queen sizes)
- No pull motors for added safety.
- Wheels - 8x heavy duty swivel castors with one lockable castor on each corner.
- Head lift angle 62° and knee break angle 42°
- Max height - 66cm from floor to top of mattress platform.
- Min height - 22cm from floor to top of mattress platform.

ACCESSORIES AVAILABLE

- High Side Rail
- Low Side Rail
- Bed Stick
- U-Assist Rail
- Over Bed Pole
- Full Length Side Rails
- Headboard
- Footboard



Available in one piece Double and Queen.