

Icare Adjustable Beds



- Adjustable beds designed for the home.
- Bed surrounds in a range of colours and fabrics.
- Optional headboards and footboards to complete the look.
- Adjustable backrest, leg rest, high / low and tilt options.
- Sizes from single through to queen.
- Partner / companion beds to match.
- Mattresses to suit all levels of pressure care needs.
- Range of side rails to assist getting in and out.

Icare Adjustable and Companion Beds

Adjustable Bed Bases

Bed	Size	Back Adjust	Leg Raise	High / Low	Tilt	Rails can be Fitted	Companion Cutout	SWL
IC111	Long Single	✓	✓	-	-	✓	✓	200kg
	King Single	✓	✓	-	-	✓	✓	200kg
	Queen – Split	✓	✓	-	-	✓	-	200kg
	King – Split	✓	✓	-	-	✓	-	200kg
IC280	Long Single	✓	✓	✓	-	✓	✓	180kg
	King Single	✓	✓	✓	-	✓	✓	180kg
IC333	Long Single	✓	✓	✓	✓	✓	✓	200kg
	King Single	✓	✓	✓	✓	✓	✓	200kg
	Double	✓	✓	✓	✓	✓	-	300kg*
	Queen	✓	✓	✓	✓	✓	-	300kg*
	Queen - Split	✓	✓	✓	✓	✓	-	300kg*
IC555	King Single	✓	✓	✓	✓	✓	✓	350kg
IC777	Long Single	✓	✓	✓	✓	✓	-	200kg
	King Single	✓	✓	✓	✓	✓	-	200kg
	Double	✓	✓	✓	✓	✓	-	300kg*
	Queen	✓	✓	✓	✓	✓	-	300kg*

* SWL of 300kg for Double, Queen and Split Queen beds is for two occupants combined, not a single occupant.

Most beds available in choice of four fabric colours and two vinyl colours (extra cost for vinyl).

IC777 and IC700 beds are Premium Homecare beds available in 3 different headboard/foot board designs and a choice of three fabric colours and three vinyl colours.

Companion (fixed) Bed Bases

Bed	Size	Back Adjust	Leg Raise	High / Low	Tilt	Rails can be Fitted	Companion Cutout	SWL
IC100	Long Single	-	-	-	-	✓	✓	200kg
	King Single	-	-	-	-	✓	✓	200kg
IC700	Long Single	-	-	-	-	✓	-	200kg
	King Single	-	-	-	-	✓	-	200kg
CB	Long Single	-	-	-	-	-	✓	100kg
	King Single	-	-	-	-	-	✓	100kg

Due to the size and weight of these beds, they cannot be installed in all locations.

Prior to order acceptance, an assessment (free of charge) must be conducted to confirm suitability.