


# Invacare Softform Mattresses

## Softform Premier:

- Deep contoured insert pad for enhanced pressure care area.
- Castellated top surface for pressure care and firm edges for stability when standing up/sitting on bed.
- Suitable for High Risk / Very High Risk when used with an adjustable bed.
- Sizes: Single and King Single.
- SWL 247kg.



## Softform Maxiglide:

- Glide system allows mattress to fully conform to the bed frame.
- Top surface remains under user throughout profiling of an electric bed.
- Castellated top surface for pressure care and firm edges for stability when standing up / sitting on bed.
- Suitable for Very High Risk when used with an adjustable bed.
- Sizes: Single and King Single.
- SWL 247kg.

