

Quick Reference Guide

Mobile Shower Commode Range

Standard tilt in space mobile shower commode (KA11507)



Standard self-propel mobile shower commode (KA120SR)



Standard transit mobile shower commode (KA114S)



Knowledge, Experience and Specialist Consultation

Years of knowledge, experience, and consultation with specialists and consumers, goes into the design and construction of K Care's premium range of shower commodes to deliver superior;

Support

Comfort

Safety

Functionality

The following quick reference guide (see over) is designed to highlight standard features, key differences and various options.

Attendant Propelled Shower Commodes

KA114SC	AP	320mm	120kg	Strap Back, Paediatric , Swing-up Arms, Swing-away Foot Rests , 4 Locking Castors
KA114SJ	AP	400mm	200kg	Strap Back, Junior with Swing-up Arms, Swing-away Foot Rests , 4 Locking Castors
KA112S	AP	445mm	200kg	Strap Back, Swing-up Arms, Sliding Foot Plate , 4 Locking Castors
KA113S	AP	445mm	200kg	Strap Back, Weight Bearing Foot Platform , 4 Locking Castors
KA114S	AP	445mm	200kg	Strap Back, Swing-up Arms, Swing-away Foot Rests , 4 Locking Castors
KA114S-HA	AP	445mm	200kg	Height Adjustable , Strap Back, Swing-up Arms, Swing-away Foot Rests , 4 Locking Castors
KA116S	AP	445mm	200kg	Strap Back, Folding with Swing-up Arms, Swing-away Foot Rests , 4 Locking Castors
KA112S50	AP	500mm	200kg	Strap Back Rest, Swing-up Padded Arms, Sliding Foot Plate , 4 Locking Castors
KA113S50	AP	500mm	200kg	Strap Back Rest, Swing-up Padded Arms, Weight Bearing Foot Platform , 4 Locking Castors
KA114S50	AP	500mm	200kg	Strap Back Rest, Swing-up Padded Arms, Swing-away Foot Rests , 4 Locking Castors
KA112S55	AP	550mm	300kg	Padded Back Rest, Removable Padded Arms, Sliding Foot Plate , 4 Locking Castors
KA113S55	AP	550mm	300kg	Padded Back Rest, Removable Padded Arms, Weight Bearing Foot Platform , 4 Locking Castors
KA114S55	AP	550mm	300kg	Padded Back Rest, Removable Padded Arms, Swing-away Foot Rests , 4 Locking Castors
KA112S60	AP	600mm	300kg	Padded Back Rest, Removable Padded Arms, Sliding Foot Plate , 4 Locking Castors
KA113S60	AP	600mm	300kg	Padded Back Rest, Removable Padded Arms, Weight Bearing Foot Platform , 4 Locking Castors
KA114S60	AP	600mm	300kg	Padded Back Rest, Removable Padded Arms, Swing-away Foot Rests , 4 Locking Castors
KA112S65	AP	650mm	300kg	Padded Back Rest, Removable Padded Arms, Sliding Foot Plate , 4 Locking Castors
KA113S65	AP	650mm	300kg	Padded Back Rest, Removable Padded Arms, Weight Bearing Foot Platform , 4 Locking Castors
KA114S65	AP	650mm	300kg	Padded Back Rest, Removable Padded Arms, Swing-away Foot Rests , 4 Locking Castors

Self Propelled Shower Commodes

KA120SRC	SPR	320mm	120kg	Strap Back, Paediatric , Swing-up Arms, Swing-away Foot Rests , 2 Locking Castors
KA120SRJ	SPR	400mm	200kg	Strap Back, Junior , Swing-up Arms, Swing-away Foot Rests , 2 Locking Castors
KA120SR	SPR	445mm	200kg	Strap Back, Swing-up Arms, Swing-away Foot Rests , 2 Locking Castors
KA120SR-HA	SPR	445mm	200kg	Height Adjustable , Strap Back, Swing-up Arms, Swing-away Foot Rests , 2 Locking Castors
KA120SF	SPF	445mm	200kg	Strap Back, Swing-up Arms, Flip-up Foot Plates , 2 Locking Castors
KA123S	SPR	445mm	200kg	Strap Back, Swing-up Arms, Weight Bearing Foot Platform , 2 Locking Castors
KA130S	SPR	445mm	200kg	Strap Back, Folding , Swing-up Arms, Swing-away Foot Rests , 2 Locking Castors
KA120SR50	SPR	500mm	200kg	Strap Back, Swing-up Padded Arms, Swing-away Foot Rests , 2 Locking Castors
KA123S50	SPR	500mm	200kg	Strap Back, Swing-up Arms, Weight Bearing Foot Platform , 2 Locking Castors
KA120SR55	SPR	550mm	300kg	Padded Back Rest, Removable Padded Arms, Swing-away Foot Rests , 2 Locking Castors
KA120SR60	SPR	600mm	300kg	Padded Back Rest, Removable Padded Arms, Swing-away Foot Rests , 2 Locking Castors
KA120SR65	SPR	650mm	300kg	Padded Back Rest, Removable Padded Arms, Swing-away Foot Rests , 2 Locking Castors

Tilt in Space Shower Commodes

KA115S40	AP	400mm	200kg	Gas Lift Tilt-in-Space, Junior , Strap Back, Swing-up Padded Arms , Padded Calf Rests, Swing-away Foot Rests, 4 Locking Castors
KA11507	AP	445mm	200kg	Gas Lift Tilt-in-Space, Strap Back, Swing-up Padded Arms , Padded Calf Rests, Swing-away Foot Rests, 4 Locking Castors
KA11507-HA	AP	445mm	200kg	Gas Lift Tilt-in-Space, Height Adjustable , Strap Back, Swing-up Padded Arms , Padded Calf Rests, Swing-away Foot Rests, 4 Locking Castors
KA11508	AP	445mm	200kg	Gas Lift Tilt-in-Space, Strap Back, Drop-in Padded Arms , Padded Side Panels , Padded Calf Rests, Swing-away Foot Rests, 4 Locking Castors
KA11512	AP	500mm	200kg	Gas Lift Tilt-in-Space, Strap Back Rest, Swing-up Padded Arms , Padded Calf Rests, Swing-away Foot Rests, 4 Locking Castors
KA11506	AP	600mm	300kg	Electric Lift Tilt-in-Space, Padded Back Rest, Removable Padded Arms, Padded Calf Rests, Swing-away Foot Rests, 4 Locking Twin Castors
KA115065	AP	650mm	300kg	Electric Lift Tilt-in-Space, Padded Back Rest, Removable Padded Arms, Padded Calf Rests, Swing-away Foot Rests, 4 Locking Twin Castors

REHAB Shower Commodes*

KA114SC-REHAB	AP	320mm	120kg	Height adj, Tension adj Back, Paediatric , Height adj Swing-up Arms, Swing-away Foot Rests , 4 Locking Castors
KA114SK-REHAB	AP	360mm	120kg	Height adj, Tension adj Back, Child , Height adj Swing-up Arms, Swing-away Foot Rests , 4 Locking Castors
KA114SJ-REHAB	AP	400mm	200kg	Height adj, Tension adj Back, Junior , Height adj Swing-up Arms, Swing-away Foot Rests , 4 Locking Castors
KA114S-REHAB	AP	445mm	200kg	Height adj, Tension adj Back, Height adj Swing-up Arms, Swing-away Foot Rests , 4 Locking Castors
KA114S50-REHAB	AP	500mm	200kg	Height adj, Tension adj Back Rest, Height adj Swing-up Arms, Swing-away Foot Rests , 4 Locking Castors
KA120SRC-REHAB	SPR	320mm	120kg	Height adj, Tension adj Back, Paediatric , Height adj Swing-up Arms and Swing-away Foot Rests , 2 Locking Castors
KA120SRK-REHAB	SPR	360mm	120kg	Height adj, Tension Adj Back, Child , Height adj Swing-up Arms and Swing-away Foot Rests , 2 Locking Castors
KA120SRJ-REHAB	SPR	400mm	200kg	Height adj, Tension adj Back, Junior , Height adj Swing-up Arms and Swing-away Foot Rests , 2 Locking Castors
KA120SR-REHAB	SPR	445mm	200kg	Height adj, Tension Adj Back, Height adj Swing-up Arms and Swing-away Foot Rests , 2 Locking Castors
KA120SR50-REHAB	SPR	500mm	200kg	Height adj, Tension Adj Back, Height adj Swing-up Padded Arms and Swing-away Foot Rests , 2 Locking Castors
KA115S40-REHAB	AP	400mm	200kg	Height adj, Gas Lift Tilt-in-Space, Tension Adj Back, Height adj Swing-up Arms, Padded Calf Rests and Swing-away Foot Rests , 4 Locking Castors
KA11507-REHAB	AP	445mm	200kg	Height adj, Gas Lift Tilt-in-Space, Tension adj Back, Height adj Swing-up Arms , Padded Calf Rests, Swing-away Foot Rests, 4 Locking Castors
KA11512-REHAB	AP	500mm	200kg	Height adj, Gas Lift Tilt-in-Space, Tension adj Back, Height adj Swing-up Arms , Padded Calf Rests, Swing-away Foot Rests, 4 Locking Castors

*All REHAB models have height adjustable frames, height adjustable arms. Plus mounting plates for – lateral thoracic pads, hip & knee guides

For more information on individual models, go to www.kcare.com.au and enter the product code into the search box

For specific customising enquiries: customcommodes@kcare.com.au

Please contact customerservice@kcare.com.au to find out more information on specific models in each range.

www.kcare.com.au

