

Height Adjustable Tilt in Space Shower Commode



30° of tilt using twin gas strut for smooth and easy tilt



Standard left and right side opening



Removable, swing-away and height adjustable footrests

Features

- Height adjustable
- Standard left and right side opening
- Extra high PVC strap backrest for improved comfort and safety
- Integrated push handle for maneuverability
- Swing-up armrest for unrestricted transfers
- Padded armrest and calf supports for improved comfort and support
- Removable upper backrest section
- Removable, swing-away and height adjustable footrests
- 30° of tilt using twin gas strut for smooth and easy tilt
- Conforms to AS/NZS 3973:2009 for safety and peace of mind

SPECIFICATIONS

Description	Max User Weight	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat frame Height Range
Height Adjustable Tilt in Space Shower Commode (KA11507-HA)	200kg	445mm	435mm	855mm	580mm	1070mm	540-465mm

Please contact customerservice@kcare.com.au for more information.

www.kcare.com.au