



# Manual Handling / Slings

Hoist/Personal Lifting Device

Sling for Hoist

Accessories



# Manual Handling / Slings

Hoist / Personal Lifting Device



## Multi Electric Lift

- With the touch of a button, is capable of lifting patients in either a standing or sitting position
- Powered lifting arm and leg gating control
- Multi-function handset make for easier, faster and safer patient lifting and transport
- Intelligent electronics within the lifter record usage data, battery charge and the condition of the battery
- Electronic overload feature prevents the lifter operating if it detects a weight above its capacity
- A maintenance required indicator prompts users to service the Multi Lift to ensure ongoing safety and extend the operating life
- The Multi Lift has many features and is compatible with a range of slings



		MAX USER WEIGHT (KG)	OVERALL LENGTH (MM)	OVERALL HEIGHT (MM)	PRODUCT WEIGHT (KG)	TURNING RADIUS (MM)	MIN LIFT HEIGHT (MM)	MAX LIFT HEIGHT (MM)	OVERALL WIDTH (MM)
Multi Electric Lift	KH401	180	1310	1220	58	650	430	1700	700



## Patient Lifter - Quick Stand

- Independently tested to conform to ISO 10535 for safety and peace of mind.
- Linak Actuation.
- Electric lifting arm and leg adjustment using handset.
- Plug in rechargeable battery system.
- Hand operated emergency lowering.
- Emergency stop button.
- Wide lateral push handles for manoeuvrability.
- Adjustable height knee pad.
- Fold up foot plate.



		MAX USER WEIGHT (KG)	OVERALL LENGTH (MM)	OVERALL HEIGHT (MM)	PRODUCT WEIGHT (KG)	MIN LIFT HEIGHT (MM)	MAX LIFT HEIGHT (MM)	OVERALL WIDTH (MM)
Patient Lifter - Quick Stand	KH400G	200	1050	910	45	820	1620	670



### Patient Lifter - Ezy Lift

- Independently tested to conform to ISO 10535 for safety and peace of mind.
- Linak Actuation.
- Electric lifting arm and leg adjustment using handset.
- Plug in rechargeable battery system.
- Hand operated emergency lowering.
- Emergency stop button.



		MAX USER WEIGHT (Kg)	OVERALL LENGTH (MM)	OVERALL HEIGHT (MM)	PRODUCT WEIGHT (Kg)	OVERALL WIDTH (MM)
Patient Lifter - Ezy Lift	KH405F	225	1300	1315	53	700

## Manual Handling / Slings

### Sling for Hoist



### Standing Transfer Sling Hook and Loop

- Ideal for transferring patients from room to room
- Polyester material
- High density foam offers excellent patient comfort and support
- Quick and easy to fit
- Neoprene waist strap with Hook and Loop
- 4 colour-coded attachment loops (Red, Yellow, Green and Blue) for different size patients
- Machine washable



		MAX USER WEIGHT (Kg)
Standing Transfer Sling, Hook and Loop Waist Strap, Medium	KH400SF	250
Standing Transfer Sling, Hook and Loop Waist Strap, Large	KH400SFL	250
Standing Transfer Sling, Hook and Loop Waist Strap, Small	KH400SFS	250



### Standing Transfer Sling, Hook and Loop and Buckle

- Ideal for transferring patients from room to room
- Polyester material
- High density foam offers excellent patient comfort and support
- Quick and easy to fit
- Neoprene waist strap with Hook and Loop
- Webbing waist strap with Buckle
- 4 colour-coded attachment loops (Red, Yellow, Green and Blue) for different size patients
- Machine washable



MAX USER WEIGHT (KG)

Standing Transfer Sling, Hook and Loop and Buckle Waist Strap, Medium	KH400SFB	250
Standing Transfer Sling, Hook and Loop and Buckle Waist Strap, Large	KH400SFLB	250



### General Purpose Sling, Curved Top

- Ideal for toileting patients
- Polyester material
- Excellent back and leg support
- Padding throughout the sling offers excellent patient comfort
- Quick and easy to fit
- Neoprene waist strap with Velcro
- 4 colour-coded attachment loops (Red, Yellow, Green and Blue) at shoulders and at legs for different size patients
- 2 handles to assist in manoeuvring patient
- Machine washable



MAX USER WEIGHT (KG)

General Purpose Sling, Curved Top with Hook and Loop Waist Strap, Poly, Medium	KH406N2	250
General Purpose Sling, Curved Top with Hook and Loop Waist Strap, Poly, Large	KH406N3	250
General Purpose Sling, Curved Top with Hook and Loop Waist Strap, Poly, Extra Large	KH406N4	250
General Purpose Sling, Curved Top with Hook and Loop Waist Strap, Poly, Small	KH406N1	250



### General Purpose Sling, Straight Top, Mesh

- Ideal for bathing patients
- Mesh/netting material
- Excellent back and leg support
- Padding in the leg area offers excellent patient comfort
- Quick and easy to fit
- 4 colour-coded attachment loops (Red, Yellow, Green and Blue) at shoulders and at legs for different size patients
- 2 handles to assist in manoeuvring patient
- Machine washable



MAX USER WEIGHT (Kg)

General Purpose Sling, Straight Top, Mesh, Small	KH406NE1	250
General Purpose Sling, Straight Top, Mesh, Medium	KH406NE2	250
General Purpose Sling, Straight Top, Mesh, Large	KH406NE3	250
General Purpose Sling, Straight Top, Mesh, Extra Large	KH406NE4	250



### General Purpose Sling, Straight Top, Poly

- Ideal for toileting patients
- Polyester material
- Excellent back and leg support
- Padding throughout the sling offers excellent patient comfort
- Quick and easy to fit
- 4 colour-coded attachment loops (Red, Yellow, Green and Blue) at shoulders and at legs for different size patients
- 2 handles to assist in manoeuvring patient
- Machine washable



MAX USER WEIGHT (Kg)

General Purpose Sling, Straight Top, Poly, Small	KH406NF1	250
General Purpose Sling, Straight Top, Poly, Medium	KH406NF2	250
General Purpose Sling, Straight Top, Poly, Large	KH406NF3	250
General Purpose Sling, Straight Top, Poly, Extra Large	KH406NF4	250



### General Purpose Sling with Head Support, Mesh

- Ideal for bathing patients
- Mesh/netting material
- Excellent back, head and leg support, complete with 6 flexible strengthening inserts in the back area
- Padding in the back and leg areas offer excellent patient comfort
- Quick and easy to fit
- 4 colour-coded attachment loops (Red, Yellow, Green and Blue) at shoulders and at legs for different size patients
- 3 handles to assist in manoeuvring patient
- Machine washable



MAX USER WEIGHT (KG)

General Purpose Sling with Head Support, Mesh, Small	KH407NE1	250
General Purpose Sling with Head Support, Mesh, Medium	KH407NE2	250
General Purpose Sling with Head Support, Mesh, Large	KH407NE3	250
General Purpose Sling with Head Support, Mesh, Extra Large	KH407NE4	250



### Full Body Sling with Divided Leg, Poly

- Ideal for toileting patients that require full body support
- Polyester material
- Divided leg support for greater support and to create a commode aperture
- Full body and head support, complete with 2 flexible strengthening inserts in the back area
- Padding throughout offers excellent patient comfort
- Quick and easy to fit
- 4 colour-coded attachment loops (Red, Yellow, Green and Blue) at shoulders and at legs for different size patients
- 3 handles to assist in manoeuvring patient
- Machine washable



MAX USER WEIGHT (KG)

Full Body Sling with Divided Leg, Poly, Small	KH408N1	250
Full Body Sling with Divided Leg, Poly, Medium	KH408N2	250
Full Body Sling with Divided Leg, Poly, Large	KH408N3	250
Full Body Sling with Divided Leg, Poly, Extra Large	KH408N4	250




### Full Body Sling with Divided Leg, Mesh

- Ideal for bathing patients that require full body support
- Mesh/netting material
- Divided leg support for greater support and to create a commode aperture
- Full body and head support, complete with 2 flexible strengthening inserts in the back area
- Padding throughout offers excellent patient comfort
- Quick and easy to fit
- 4 colour-coded attachment loops (Red, Yellow, Green and Blue) at shoulders and at legs for different size patients
- 3 handles to assist in manoeuvring patient
- Machine washable



MAX USER WEIGHT (KG)

Full Body Sling with Divided Leg, Mesh, Small	KH408NE1	250
Full Body Sling with Divided Leg, Mesh, Medium	KH408NE2	250
Full Body Sling with Divided Leg, Mesh, Large	KH408NE3	250
Full Body Sling with Divided Leg, Mesh, Extra Large	KH408NE4	250




### Hammock Sling, Poly

- Ideal for toileting patients that require full body support
- Polyester material
- Commode aperture
- Full body and head support, complete with 2 large flexible strengthening inserts in the back area
- Padding in the back area offers excellent patient comfort
- Quick and easy to fit
- 4 colour-coded attachment loops (Red, Yellow, Green and Blue) at shoulders and at legs for different size patients
- 2 handles to assist in manoeuvring patient
- Machine washable



MAX USER WEIGHT (KG)

Hammock Sling, Poly, Small	KH409N1	250
Hammock Sling, Poly, Medium	KH409N2	250
Hammock Sling, Poly, Large	KH409N3	250
Hammock Sling, Poly, Extra Large	KH409N4	250



### Hammock Sling, Mesh

- Ideal for bathing patients that require full body support
- Mesh/netting material
- Commode aperture
- Full body and head support, complete with 2 large flexible strengthening inserts in the back area
- Padding in the back area offers excellent patient comfort
- Quick and easy to fit
- 4 colour-coded attachment loops (Red, Yellow, Green and Blue) at shoulders and at legs for different size patients
- 2 handles to assist in manoeuvring patient
- Machine washable



		MAX USER WEIGHT (KG)
Hammock Sling, Mesh, Small	KH409NE1	250
Hammock Sling, Mesh, Medium	KH409NE2	250
Hammock Sling, Mesh, Large	KH409NE3	250



### Quick Fit General Purpose Sling, Poly

- Ideal for toileting patients that require full body support
- Polyester material
- Divided leg support
- Full body and head support
- Padding in the leg area offers excellent patient comfort
- Quick and easy to fit
- 4 colour-coded attachment loops (Red, Yellow, Green and Blue) at shoulders and at legs for different size patients
- 3 positioning handles
- Machine washable



		MAX USER WEIGHT (KG)
Quick Fit General Purpose Sling, Poly, Small	KH410N1	250
Quick Fit General Purpose Sling, Poly, Medium	KH410N2	250
Quick Fit General Purpose Sling, Poly, Large	KH410N3	250
Quick Fit General Purpose Sling, Poly, Extra Large	KH410N4	250





### General Purpose Sling with Head Support, Poly.

- Ideal for bathing patients with a mesh/netting material
- Easy to fit to lifters with 4 colour coded attachment loops
- A 10mm thick padding in the back, leg and head area for excellent patient comfort
- Excellent back and head support with flexible strengthening inserts in the back area
- Commode aperture
- 3 positioning straps/handles to assist with manoeuvring patient
- Machine washable



MAX USER WEIGHT (KG)

General Purpose Sling with Head Support, Poly, Small	KH407N1	250
General Purpose Sling with Head Support, Poly, Medium	KH407N2	250
General Purpose Sling with Head Support, Poly, Large	KH407N3	250
General Purpose Sling with Head Support, Poly, Extra Large	KH407N4	250

## Manual Handling / Slings

Accessories



### Pivot Frame

- Pivot Frame with hook or key hole location points for Kerry sitting patient lifters
- Offers easy patient positioning
- Easy pivot action allows patients to be placed into a lying or seating position



Pivot Frame with Hook Location Points for Kerry Patient Lifters	KH40PFH
Pivot Frame with Key Hole Location Points for Kerry Sitting Patient Lifters	KH40PF



### Digital Lift Scale for Kerry Patient Lifters

- Digital Lift Scale for Kerry Patient Lifters
- Easily installed between boom and yoke
- Weigh patients up to 400KG



Digital Lift Scale for Kerry Patient Lifters	KH40EWS
--	---------