

# Mobile Shower Commode | KA114S

## Push Chair - Attendant Propelled

### 445mm Wide Seat



The KA114S Mobile Shower Commode – Push Chair Attendant propelled is designed for multi-user situations found in hospitals, clinics and nursing homes.

#### FEATURES

Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.

3 strap back rest for comfort and easy access.

Integrated push-handle for manoeuvrability.

Swing-up arms for unrestricted side transfer.

Removable swing-away adjustable height foot rests for improved comfort and easy access.

Large castors make the chair easy to manoeuvre.

Durable and long lasting with a stainless steel frame, impervious PVC straps and water resistant castors.

Designed to be used with a K Care patented removable seat (optional KA11CPU seat shown).

Please contact K Care for information on product options available and any customisation requirements.

#### SPECIFICATIONS

Overall Width: **520mm**  
Overall Depth: **930mm**  
Overall Height: **985mm**

Seat Width: **445mm**  
Seat Depth: **430mm**  
Seat Height: **590mm** (floor to top of 40mm seat)

Backrest Height: **420mm** (from top of 40mm seat)

Castors: **125mm Total Lock (Front)**  
**125mm Swivel (Rear)**

#### Materials

Frame: **Stainless Steel**  
Frame Finish: **Electro Polished**  
Straps: **PVC Plastic**  
Castors: **Plastic (Tente)**

Product Weight: **10kg** (approx.)

Max. User Weight: **160kg**

ARTG No.: **209826**

**Please note:**  
Specifications are subject to change without notice.  
Dimensions shown are nominal.

JAN 2017