

SELF PROPELLED

Height Adjustable Shower Commode

KA120SR-HA

FEATURES

- 445mm seat width
- Height adjustability of the frame
- Flip up arms for unrestricted side transfer
- Strap back rest for comfort and easy access
- Durable and long lasting with a stainless steel frame
- Swing away footrests
- Large self-propelling rear wheels are flush with seat top for unrestricted side transfer
- 2 x Front 125mm locking swivel castors
- Conforms to AS/NZS 3973:2009 for safety and peace of mind



SPECIFICATIONS

Max User Weight	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat Frame Height Front Range	Floor to Seat Frame Height Rear Range
200kg	445mm	430mm	400mm	680mm	1050mm	550 - 475mm	527 - 452mm

Improving your
life everyday™

kcare.com.au

KCare
Healthcare Solutions

