

Mobile Shower Commode Push Chair - Self Propelled 445mm Wide Seat

KA120SR



The KA120SR Mobile Shower Commode – Push Chair Self propelled is designed for multi-user situations found in hospitals, clinics and nursing homes.

FEATURES

Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.

3 strap back rest for comfort and easy access.

Swing-up arms for unrestricted side transfer.

Large self-propelling rear wheels are flush with seat top for unrestricted side transfer.

Removable swing-away adjustable height foot rests for improved comfort and easy access.

Three wheel positions for flexible set-up.

Durable and long lasting with a stainless steel frame, impervious PVC straps and water resistant castors.

Designed to be used with a K Care patented removable seat (optional KA11CPU seat shown).

Please contact K Care for information on product options available and any customisation requirements.

SPECIFICATIONS

Overall Width: **615mm**
Overall Depth: **1050mm**
Overall Height: **985mm**

Seat Width: **445mm**
Seat Depth: **430mm**
Seat Height: **590mm** (floor to top of 40mm seat)

Backrest Height: **420mm** (from top of 40mm seat)

Wheels: **560mm (Rear)**
Castors: **125mm Swivel (Front)**

Materials
Frame: **Stainless Steel**
Frame Finish: **Electro Polished**
Straps: **PVC Plastic**
Castors: **Plastic (Tente)**

Product Weight: **15kg** (approx.)

Max. User Weight: **160kg**

ARTG No.: **209826**

Please note:
Specifications are subject to change without notice.
Dimensions shown are nominal.

JAN 2017