

Karma Flexx Wheelchair



Adjustable wheelchair:

- Folding light weight aluminium wheelchair.
- Swing away, removable leg rests. Swing in / out.
- Flip-up, height adjustable armrests with ergonomic pads.
- Quick release wheels – 24” – with solid tyres.
- Solid 7” castors – available in different sizes.
- Push to lock brakes.
- Adjustable seat height.
- Adjustable seat depth.
- Adjustable backrest angle.
- Adjustable centre of gravity.
- Extensive range of options and accessories.
- Warranty - 10 years on frame.
- SWL 130kg – standard.
- SWL 170kg – heavy duty.

Size	Seat Width	Seat Depth	Seat Height	Back Height	SWL	Chair Weight
16” x 16”	41cm	41cm	50cm	43cm	130kg	14.4kg
18” x 16”	46cm	41cm	50cm	43cm	130kg	14.6kg
18” x 18”	46cm	46cm	50cm	43cm	130kg	14.9kg
18” x 20”	46cm	52cm	50cm	43cm	130kg	15.2kg
20” x 18” HD	51cm	46cm	52cm	43cm	170kg	17kg
22” x 18” HD	55cm	46cm	52cm	43cm	170kg	18kg
24” x 18” HD	60.5cm	46cm	52cm	43cm	170kg	19kg

Available Accessories

Elevating Leg Rest (each)	Seat Extension Kit
Stump Support (each)	Tie-Down Bracket Kit
Vertical Leg Rests (each)	IV Pole
Angle Adjustable Footplate (each)	Trough Arm
Anti-tip Bars and Wheels (pair)	Oxygen Bottle Holder
Tension Adjustable Backrest	Height Adjustable Push Handles (pair)
Attendant Brakes (pair)	Height Adjustable Backrest Upholstery and Adjustable Backrest Canes
Folding Headrest	Backpack
Armrest Pouch Bag	Rear Carry Bag