

Karma VIP2 Wheelchair

Adjustable wheelchair:

- Folding tilt-in-space wheelchair.
- Self-propel wheels – 20in diameter.
- Seat tilt from 0 to 30 degrees.
- Backrest angle from 0 to 35 degrees.
- Tension-adjustable backrest.
- Curved headrest, removable.
- Flip-up, height adjustable armrests with ergonomic pads.
- Swing away, removable leg rests.
- Optional elevating leg rests.
- Optional vertical leg rests.
- Optional stump supports.
- Size when folded: 84cm x 70cm x 42cm.



Folds for transporting

Size: 84cm wide x 70cm high

(note: Transit wheel model shown not available in Australia)

Size	Seat Width	Seat Depth	Seat Height	Back Height	SWL	Chair Weight
16" x 16"	40cm	41-46cm	49cm	49cm	115kg	23.4kg
18" x 16"	45cm	41-46cm	49cm	49cm	115kg	23.4kg

Available Accessories

- Elevating Leg Rests
- Stump Support
- Vertical Leg Rests
- Trough Arms