




Transit	Karma Ergo Lite	Karma Ergo Lite Deluxe	Karma S-Ergo125																														
Image:																																	
Description:	<p>One of the lightest transit wheelchairs available, easily folded for travelling, ideal for regular short trips to the shops or doctor. Note: Leg rests are fixed – do not swing away and cannot be removed.</p>	<p>One of the lightest transit wheelchairs available, easily folded for travelling, ideal for regular short trips to the shops or doctor. Comes with swing away, removable leg rests.</p>	<p>A lightweight folding wheelchair with an ergonomic seat that helps with pressure care and prevents the user from sliding forward when sitting, suitable for all-day use.</p>																														
Key Features:	<ul style="list-style-type: none"> • Ultra-lightweight aluminium frame • Footplates are height adjustable • Padded arm rests • Ergonomic padded seat and back • Fold down back rest • Attendant hand brakes • Easy folding for storage and transportation 	<ul style="list-style-type: none"> • Ultra-lightweight aluminium frame • Footplates are detachable and height adjustable • Padded arm rests • Ergonomic padded seat and back • Fold down back rest • Attendant hand brakes • Easy folding for storage and transportation 	<ul style="list-style-type: none"> • Lightweight aluminium frame • Footplates are removable and height-adjustable • Padded swing-back armrests • Ergonomic padded seat and back • Fold down back rest • Attendant hand brakes and manual brakes • Easy folding for storage and transportation 																														
Models:	<table border="1"> <thead> <tr> <th><u>W x D</u></th> <th><u>Net wt.</u></th> <th><u>SWL</u></th> </tr> </thead> <tbody> <tr> <td>40cm x 40cm</td> <td>8.1kg</td> <td>100kg</td> </tr> <tr> <td>46cm x 40cm</td> <td>8.3kg</td> <td>100kg</td> </tr> </tbody> </table>	<u>W x D</u>	<u>Net wt.</u>	<u>SWL</u>	40cm x 40cm	8.1kg	100kg	46cm x 40cm	8.3kg	100kg	<table border="1"> <thead> <tr> <th><u>W x D</u></th> <th><u>Net wt.</u></th> <th><u>SWL</u></th> </tr> </thead> <tbody> <tr> <td>40cm x 40cm</td> <td>8.7kg</td> <td>100kg</td> </tr> <tr> <td>46cm x 40cm</td> <td>8.8kg</td> <td>100kg</td> </tr> </tbody> </table>	<u>W x D</u>	<u>Net wt.</u>	<u>SWL</u>	40cm x 40cm	8.7kg	100kg	46cm x 40cm	8.8kg	100kg	<table border="1"> <thead> <tr> <th><u>W x D</u></th> <th><u>Net wt.</u></th> <th><u>SWL</u></th> </tr> </thead> <tbody> <tr> <td>40cm x 43cm</td> <td>13.1kg</td> <td>115kg</td> </tr> <tr> <td>46cm x 43cm</td> <td>13.1kg</td> <td>115kg</td> </tr> <tr> <td>51cm x 43cm</td> <td>13.3kg</td> <td>115kg</td> </tr> </tbody> </table>	<u>W x D</u>	<u>Net wt.</u>	<u>SWL</u>	40cm x 43cm	13.1kg	115kg	46cm x 43cm	13.1kg	115kg	51cm x 43cm	13.3kg	115kg
<u>W x D</u>	<u>Net wt.</u>	<u>SWL</u>																															
40cm x 40cm	8.1kg	100kg																															
46cm x 40cm	8.3kg	100kg																															
<u>W x D</u>	<u>Net wt.</u>	<u>SWL</u>																															
40cm x 40cm	8.7kg	100kg																															
46cm x 40cm	8.8kg	100kg																															
<u>W x D</u>	<u>Net wt.</u>	<u>SWL</u>																															
40cm x 43cm	13.1kg	115kg																															
46cm x 43cm	13.1kg	115kg																															
51cm x 43cm	13.3kg	115kg																															
Popular accessories:	<p>Folding headrest Backpack Rear carry bag Armrest pouch bag Anti-tip bars pair Taller Push Handles pair (7cm higher)</p>	<p>Folding headrest Backpack Rear carry bag Armrest pouch bag Anti-tip bars pair Taller Push Handles pair (7cm higher)</p>	<p>Folding Headrest Backpack Rear carry bag Armrest pouch bag Stump support Elevating leg rest Right / Left Anti-tip bars pair Taller Push Handles pair (7cm higher)</p>																														

Self-propel	Karma Ergo Lite Deluxe	Karma S-Ergo 125	Karma Flexx																																										
Image:																																													
Description:	One of the lightest self-propelling wheelchairs available, easily folded for travelling, ideal for frequent short trips to the shops or doctor.	A lightweight folding wheelchair with an ergonomic seat that helps with pressure care and prevents the user from sliding forward when sitting, suitable for all-day use.	A configurable self-propel wheelchair, suitable for all-day use as it can be adjusted in several different ways to accomplish optimum comfort and performance for users.																																										
Key Features:	<ul style="list-style-type: none"> • Ultra-lightweight aluminium frame • Footplates are detachable and height adjustable • Padded arm rests • Ergonomic padded seat and back • Fold down back rest • Attendant brakes and manual brakes • Easy folding for storage 	<ul style="list-style-type: none"> • Lightweight aluminium frame • Footplates are removable and height-adjustable • Padded swing-back armrests • Ergonomic padded seat and back • Fold down back rest • Quick release wheels with solid rear tyres • Larger rear wheels are easier to reach and make self-propelling easier • Attendant brakes and manual brakes 	<ul style="list-style-type: none"> • Lightweight aluminium frame • Footplates are removable and height-adjustable • Padded armrests are removable and height-adjustable • Padded backrest is angle adjustable • Padded seat is height adjustable and a range of accessories can adjust seat depth • Manual brakes 																																										
Models:	<table border="1"> <thead> <tr> <th>W x D</th> <th>Net wt.</th> <th>SWL</th> </tr> </thead> <tbody> <tr> <td>40cm x 40cm</td> <td>10.2kg</td> <td>100kg</td> </tr> <tr> <td>46cm x 40cm</td> <td>10.3kg</td> <td>100kg</td> </tr> </tbody> </table>	W x D	Net wt.	SWL	40cm x 40cm	10.2kg	100kg	46cm x 40cm	10.3kg	100kg	<table border="1"> <thead> <tr> <th>W x D</th> <th>Net wt.</th> <th>SWL</th> </tr> </thead> <tbody> <tr> <td>40cm x 43cm</td> <td>14.2kg</td> <td>115kg</td> </tr> <tr> <td>46cm x 43cm</td> <td>14.5kg</td> <td>115kg</td> </tr> <tr> <td>51cm x 43cm</td> <td>14.8kg</td> <td>115kg</td> </tr> </tbody> </table>	W x D	Net wt.	SWL	40cm x 43cm	14.2kg	115kg	46cm x 43cm	14.5kg	115kg	51cm x 43cm	14.8kg	115kg	<table border="1"> <thead> <tr> <th>W x D</th> <th>Net wt.</th> <th>SWL</th> </tr> </thead> <tbody> <tr> <td>40cm x 40cm</td> <td>14.4kg</td> <td>130kg</td> </tr> <tr> <td>46cm x 40cm</td> <td>14.6kg</td> <td>130kg</td> </tr> <tr> <td>46cm x 46cm</td> <td>14.9kg</td> <td>130kg</td> </tr> <tr> <td>51cm x 46cm</td> <td>17kg</td> <td>170kg</td> </tr> <tr> <td>56cm x 46cm</td> <td>18kg</td> <td>170kg</td> </tr> <tr> <td>61cm x 46cm</td> <td>19kg</td> <td>170kg</td> </tr> </tbody> </table>	W x D	Net wt.	SWL	40cm x 40cm	14.4kg	130kg	46cm x 40cm	14.6kg	130kg	46cm x 46cm	14.9kg	130kg	51cm x 46cm	17kg	170kg	56cm x 46cm	18kg	170kg	61cm x 46cm	19kg	170kg
W x D	Net wt.	SWL																																											
40cm x 40cm	10.2kg	100kg																																											
46cm x 40cm	10.3kg	100kg																																											
W x D	Net wt.	SWL																																											
40cm x 43cm	14.2kg	115kg																																											
46cm x 43cm	14.5kg	115kg																																											
51cm x 43cm	14.8kg	115kg																																											
W x D	Net wt.	SWL																																											
40cm x 40cm	14.4kg	130kg																																											
46cm x 40cm	14.6kg	130kg																																											
46cm x 46cm	14.9kg	130kg																																											
51cm x 46cm	17kg	170kg																																											
56cm x 46cm	18kg	170kg																																											
61cm x 46cm	19kg	170kg																																											
Popular accessories:	<ul style="list-style-type: none"> Folding headrest Backpack Rear carry bag Armrest pouch bag Anti-tip bars pair Taller Push Handles pair (7cm higher) 	<ul style="list-style-type: none"> Folding headrest Backpack Rear carry bag Armrest pouch bag Stump support Elevating leg rest Right / Left Anti-tip bars pair Taller Push Handles pair (7cm higher) 	<ul style="list-style-type: none"> Elevating leg rest Right / Left Stump support Right / Left Tension adjustable backrest Anti-tip bars pair Height adjustable push handle pair <p>Note: Many more options / customisations available. Ask for detailed flyer.</p>																																										