

# Karma Wheelchairs

## Karma Ergo Lite Deluxe



Ultra lightweight wheelchair suitable for regular use and for regularly transporting in vehicles.

Features	<ul style="list-style-type: none"> <li>• Self-propel or Attendant-propel models available.</li> <li>• Desk arms</li> <li>• Fixed armrests</li> <li>• Fold down back rest</li> <li>• Swing away, removable leg rests</li> <li>• Attendant brakes and user brakes.</li> </ul>
Seat Sizes (width x depth)	41cm x 42cm 45cm x 42cm
Weight	8.8kg (Att. Propel) 10.3kg (Self Propel)
SWL	100kg
Accessories	<ul style="list-style-type: none"> <li>• Folding Headrest</li> <li>• Backpack</li> <li>• Rear Carry Bag</li> <li>• Armrest Pouch Bag</li> <li>• Anti-Tip Bars (pair)</li> <li>• Taller Push Handles (7cm higher)</li> </ul>

## Karma S-Ergo 125



High quality wheelchair with a range of accessories and options. Suitable for all day use.

Features	<ul style="list-style-type: none"> <li>• Self-propel or Attendant-propel models available</li> <li>• Desk arms</li> <li>• Flip up armrests</li> <li>• Fold down back rest</li> <li>• Quick-release wheels</li> <li>• Swing away, removable leg rests</li> <li>• Hand brakes, two-way brakes</li> <li>• Lever brakes</li> </ul>
Seat Sizes (width x depth)	41cm x 43cm 45cm x 43cm 51cm x 43cm
Weight	13.3kg (Att. Propel) 14.2kg – 14.8kg (Self Propel)
SWL	115kg
Accessories	<ul style="list-style-type: none"> <li>• Folding Headrest</li> <li>• Elevating Leg Rests</li> <li>• Stump Support, universal</li> <li>• Anti-Tip Bars</li> <li>• Taller Push Handles (7cm higher)</li> </ul>