


# Oska Mattresses



## Oska Fika:

- Multi-layered foam mattress
- Includes tapered Offloader heel relief section to reduce pressure
- Deep
- Castellated
- Suitable for High Risk when used with a profiling bed.
- Sizes: Single, King Single, Queen
- SWL 159Kg.
- Up to category 2 pressure care



## Oska Aria:

- Hybrid foam and pumpless air combination mattress
- Individual foam segments encased in interconnected air cells to provide excellent pressure care
- Includes tapered Offloader heel relief section to reduce pressure
- Firm edges for stability when stand/sitting on bed
- Suitable for Very High Risk when used with a profiling bed.
- Sizes: Single, King Single, Queen
- SWL 230Kg.
- Up to category 3 pressure care