

Ellipse Indoor Walker - Dark Blue

Product Code: 9346376080088

Features:

- AU STD. Tested 492449
- Made of Lightweight Aluminium
- Adjust Height Tool-Free
- Removable Tray
- Easy to Fold
- Easily Installable Nylon Bag



Length (cm)	66	Folded Width (cm)	26	Front Castor Diameter (cm)	17.5/ 7"
Width (cm)	54	Tray Depth (cm)	29	Rear Castor Diameter (cm)	17.5/ 7"
Max Height (cm)	94.5	Tray Width (cm)	39	Max Load (kg)	100
Minimum Height (cm)	83.5	Tray Height from the Ground (cm)	68	Net Weight (kg)	6.6