













# Ellipse Lite Aluminium Rollator

Product Code: 9346376000543 (Medium, Red)

## Features:

- AU STD. Tested 492485
- Adjust Height Tool-Free
- Adjust Brake Tool-Free
- Foldable Waterproof Seat
- Easily Installable Nylon Bag



Length (cm)	 59.5	Folded Width (cm)	 26	Front Castor Diameter (cm)	 24.2/9.5"
Width (cm)	 49	Seat Depth (cm)	 22	Rear Castor Diameter (cm)	 20.3/8"
Max Height (cm)	 97.8	Seat Width (cm)	 45.4	Max Load (kg)	 136
Minimum Height (cm)	 84	Seat Height from the Ground (cm)	 52.3	Net Weight (kg)	 8.2