













Ellipse Lite Aluminium Rollator

Product Code: 9346376000536 (Small, Red)

Features:

- AU STD. Tested 492485
- Adjust Height Tool-Free
- Adjust Brake Tool-Free
- Foldable Waterproof Seat
- Easily Installable Nylon Bag



Length (cm)	 59.5	Folded Width (cm)	 26	Front Castor Diameter (cm)	 24.2/9.5"
Width (cm)	 49	Seat Depth (cm)	 22	Rear Castor Diameter (cm)	 20.3/8"
Max Height (cm)	 94.5	Seat Width (cm)	 45.4	Max Load (kg)	 136
Minimum Height (cm)	 82.1	Seat Height from the Ground (cm)	 49	Net Weight (kg)	 8