

# Mighty Mack Rollator

Product Code: C4216-B-FX













## Features:

- ▶ AU STD. Tested 491883
- ▶ Adjust Height Tool-Free
- ▶ Adjust Brake Tool-Free
- ▶ Foldable
- ▶ Easily Installable Metal basket



## ARTG Listing:

- ▶ 227599

Length (cm)		67	Folded Width (cm)		29	Front Castor Diameter (cm)		20.3/8"
Width (cm)		63	Seat Depth (cm)		35	Rear Castor Diameter (cm)		20.3/8"
Max Height (cm)		91	Seat Width (cm)		51	Max Load (kg)		180
Minimum Height (cm)		N/A	Seat Height from the Ground (cm)		55	Net Weight (kg)		12.6