

# Mini Mack Rollator

Product Code: C4205C-FX













## Features:

- AU STD. Tested 491907
- Adjust Height Tool-Free
- Adjust Brake Tool-Free
- Foldable
- Easily Installable Metal basket



## ARTG Listing:

- 227599

Length (cm)		62	Folded Width (cm)		28	Front Castor Diameter (cm)		20.3/8"
Width (cm)		59	Seat Depth (cm)		32	Rear Castor Diameter (cm)		20.3/8"
Max Height (cm)		86	Seat Width (cm)		43	Max Load (kg)		225
Minimum Height (cm)		N/A	Seat Height from the Ground (cm)		50	Net Weight (kg)		9.5