

# Star 6" Adjustable Rollator Red

Product Code: V4261-FX Gtin: 9346376000857












## Features:

- ▶ Quick-Fit adjustable design
- ▶ 6 holes handle adjustment
- ▶ Easy to pack up for travel
- ▶ 3 seat height adjustment
- ▶ Compact Carton packaging
- ▶ Locking pin to keep walker folded
- ▶ Aust. STD Tested 493360



## ARTG Listing:

- ▶ 227600

Length (cm)		65.5-73.1	Carton Size (cm)	37*21*16	Front Castor Diameter (cm)		15.2/6"	
Width (cm)		61.7	Seat Depth (cm)		34.8	Rear Castor Diameter (cm)		15.2/6"
Max Height (cm)		93.5	Seat Width (cm)		37	Max Load (kg)		150
Minimum Height (cm)		81.5	Seat Height from the Ground (cm)		49.3-57	Net Weight (kg)		7.3