

Star 8" HD Adjustable Rollator Red

Product Code: V4264FX GTin: 9346376000864












Features:

- ▶ Quick-Fit adjustable design
- ▶ 6 holes handle adjustment
- ▶ Easy to pack up for travel
- ▶ 3 seat height adjustment
- ▶ Compact Carton packaging
- ▶ Locking pin to keep walker folded
- ▶ Aust. STD Tested 493360



ARTG Listing:

- ▶ 227600

Length (cm)		66.8-72	Carton Size (cm)	43*23*16	Front Castor Diameter (cm)		20/8"	
Width (cm)		70.8	Seat Depth (cm)		38	Rear Castor Diameter (cm)		20/8"
Max Height (cm)		96.8	Seat Width (cm)		46	Max Load (kg)		204
Minimum Height (cm)		80.8	Seat Height from the Ground (cm)		50.8-55.8	Net Weight (kg)		9.5