




REDGUM BRAND
Daily Living and
Lifestyle Products




RAMPS RANGE

ASK FOR THE BRAND

redgumbrand.com.au

SINGLE FOLDING ACCESS RAMPS



RAM2

Overall Dimensions: 600mm (L) x 725mm (W)

Product Weight: 4.8kg

Max. Height: 152mm

RAM3

Overall Dimensions: 900mm (L) x 725mm (W)

Product Weight: 6.8kg

Max. Height: 228mm

RAM4

Overall Dimensions: 1200mm (L) x 725mm (W)

Product Weight: 9.6kg

Max. Height: 305mm

RAM5

Overall Dimensions: 1500mm (L) x 725mm (W)

Product Weight: 10.3kg

Max. Height: 381mm

RAM6

Overall Dimensions: 1800mm (L) x 725mm (W)

Product Weight: 15.2kg

Max. Height: 458mm

RAM7

Overall Dimensions: 2100mm (L) x 725mm (W)

Product Weight: 10kg

Max. Height: 610mm

Maximum load on all ramps: **270kg**

Maximum angle on all ramps: **15°**

FEATURES


- Strong, lightweight aluminium construction
- Raised non-slip surface
- Folds in half for transit and storage
- Carry handle

**FOLDS IN HALF
FOR TRANSIT
AND STORAGE**





ONE PIECE RAMPS/RAMP PAIR



RAM32

Overall Dimensions: 600mm (L) x 725mm (W)
Product Weight: 4.1kg


RAM33

Overall Dimensions: 900mm (L) x 725mm (W)
Product Weight: 6.5kg

RAM34


Overall Dimensions: 1200mm (L) x 725mm (W)
Product Weight: 9.3kg

Maximum load on all ramps: **270kg**



RAM40

Overall Dimensions:
330mm (L) x 810mm (W)
Product Weight: 2kg
Maximum load: 270kg



RAM20B

Overall Dimensions: 1500mm (L) x 190mm (W)
Product Weight: 3kg (each)
Maximum load: 270kg

TELESCOPIC RAMPS

RAM70

Overall Dimensions: 1060 - 1800mm (L) x 190mm (W)
Product Weight: 3.5kg (each)

RAM74T

Overall Dimensions: 700 - 1200mm (L) x 190mm (W)
Product Weight: 4kg (each)

RAM76T

Overall Dimensions: 1100 - 1800mm (L) x 190mm (W)
Product Weight: 5.4kg (each)

RAM78T

Overall Dimensions: 1250 - 2400mm (L) x 190mm (W)
Product Weight: 6.1kg (each)




RAM80T - Triple Telescopic Ramp

Overall Dimensions: 1050 - 2400mm (L) x 190mm (W)
Product Weight: 7kg (each)

Maximum load on all ramps: **270kg**

FEATURES

- Ramp can be set at any length within range
- Simple to operate extension lock
- Portable, lightweight and easy to carry
- Strong aluminium construction
- Raised non-slip surface
- Safety edges to prevent rolling off
- Carry handle
- Slides in for transit and storage



RAMF

Overall Dimensions:
2150 (L) x 230mm (W)
Product Weight: 6.6kg (each)
Max. Load: 340kg (over pair)



RUBBER WEDGE RAMPS

RAMR1415

Overall Dimensions: 900mm (W) x 140mm (L)
Wedge height: 15mm
Product Weight: 1.6kg

RAMR1520

Overall Dimensions: 900mm (W) x 150mm (L)
Wedge height: 20mm
Product Weight: 2.1kg

RAMR1525

Overall Dimensions: 900mm (W) x 150mm (L)
Wedge height: 25mm
Product Weight: 2.5kg

RAMR2030

Overall Dimensions: 900mm (W) x 200mm (L)
Wedge height: 30mm
Product Weight: 4.2kg

RAMR2035

Overall Dimensions: 900mm (W) x 200mm (L)
Wedge height: 35mm
Product Weight: 4.7kg

RAMR2040

Overall Dimensions: 900mm (W) x 205mm (L)
Wedge height: 40mm
Product Weight: 4.9kg

RAMR2045

Overall Dimensions: 900mm (W) x 255mm (L)
Wedge height: 45mm
Product Weight: 7.2kg

FEATURES


- Solid black rubber construction
- Raised non-slip surface provides excellent grip
- Suitable for indoor and outdoor use
- Can be cut down or grooved to suit application



Maximum load on all ramps: **270kg**

REDGUM BRAND
Daily Living and
Lifestyle Products

**YOU DESERVE
THE BEST!**



RAMR2050

Overall Dimensions: 1105mm (W) x 410mm (L)
Wedge height: 50mm
Product Weight: 20kg

RAMR20100

Overall Dimensions: 1105mm (W) x 410mm (L)
Wedge height: 100mm
Product Weight: 25kg

Maximum load on all ramps: **270kg**



Please contact REDGUM Brand if you require further information on this product.

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand: 1 Business way Malaga Western Australia 6090 Part of the For-de Group Pty Ltd

ASK FOR THE BRAND

redgumbrand.com.au