

Redgum Crystal V2.0 Lightweight Wheelchair



Very light wheelchair – only 12.5Kg. Easy to transport as it has removable rear wheels and fold-down push handles.

- Lightweight Self-Propel Wheelchair (large rear wheels).
- Aluminium folding frame – colour: Silver
- Easy to remove drive wheels.
- Removable, swing-away leg rests.
- Padded arm rests.
- Fold down push handles (reduces space when transporting in car).
- Solid tyres – no punctures or flats.
- Manual brakes on rear wheels.
- Attendant brakes.
- Rear storage pocket on back of seat.
- Seat width: 470mm (18in)
- Seat depth: 420mm.
- Seat height: 520mm
- SWL: 127Kg.
- Wheelchair weight: 12.5Kg.

