

Revolution Chair



Bariatric Model Available.

Revolution Chair – Swivel & Slide Seat

- Legs are fixed in position, seat swivels and slides to enable access to a table without assistance.
- Fixed legs with long base for stability.
- Seat slides backward / forward – locks in 3 places.
- Seat swivels in full circle – locks in 3 places.
- Lock / release handles on side and rear of chair.
- Safe Working Load 120Kg.
- Ideal for use in own home, care facilities, etc.
- Available in either fixed or adjustable height.

Delivery available, charges apply.

Available for Veterans via Country Care Group.

Schedule Id: AC02

Catalogue Item: REVO-HD



	Fixed	Adjustable
Seat Height	480mm	450mm – 530mm
Seat Width	450mm	450mm
Seat Depth	485mm	485mm

Colours and Finishes

- Standard frame colour: Black satin.
- Seat and backrest available in 26 colours – Terra Nova Neo (call for colour samples chart).