

Revolution Chair – Bariatric

Revolution Chair – Swivel & Slide Seat



- Legs are fixed in position, seat swivels and slides to enable access to a table without assistance.
- Fixed legs with long base for stability.
- Seat slides backward / forward – locks in 3 places.
- Seat swivels in full circle.
- Lock / release handle on side of chair.
- Safe Working Load 250kg.
- Ideal for use in own home, care facilities, etc.



← Seat Slides →



↻ Seat Swivels ↻

Dimensions

Seat Height	480mm
Seat Width	660mm
Seat Depth	450mm

Colours and Finishes

- Standard frame colour: Black satin.
- Seat and backrest available in 26 colours – Terra Nova Neo (call for colour samples chart).