

# SHOPRIDER®

Mobility

## 2023 MODEL UPGRADES

RED  
9LSP

BLUE  
9LSPB

GREY  
9LSPG


# ALLROUNDER

## NEW FEATURES

- Easy access control
- User friendly tiller / dash
- USB charger
- Accessory / storage pocket
- Stylish fluted rims
- Double nudge bar
- LED indicators
- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- LED Headlight
- Rear lights
- Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket



## EASY ACCESS CONTROL



## SPECIFICATIONS

- Overall: **1300mm (L) x 640mm (W) x 1200mm (H)**
- Turning Circle: **1440mm**
- Speed: **10km/h**
- Maximum range: **Up to 32km**
- Safe climb angle: **9°**
- Maximum climb angle: **12°**
- Motor: **1.5Hp high torque motor (4 pole)**
- Batteries: **2 x 50 Amp minimum (supplied)**
- Charger: **8 Amp**
- Tyres (**Pneumatic**): **Ø330 x 100mm (4.00-5) - Front + Rear**
- Product Weight: **99kg - with batteries**
- Maximum User Weight: **159kg - on flat ground**

Please contact **SHOPRIDER** for more information on each product

**P: 08 9248 4180**  
**E: orders@for-de.com.au**  
 1 Business Way Malaga Western Australia 6090  
 Part of the **For-de Group Pty Ltd**

SHOPRIDER strives to provide the highest quality products possible that are subject to change without notice.

**10 YEAR WARRANTY**  
 terms and conditions apply

**IMPORTANT NOTE:**  
 It is the responsibility of the user to ensure the suitability of this product for its intended use.  
 The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

ASK FOR THE BRAND

shoprider.com.au