



EXPLORER


888SE

FEATURES

- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- LED Headlight
- Rear lights
- Large battery capacity
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket
- Front bumper bar

SPECIFICATIONS

Overall: 1280mm (L) x 600mm (W) x 1070mm (H)
 Turning Circle: 1500mm
 Speed: 8km/h
 Maximum range: Up to 40km
 Safe climb angle: 9°
 Maximum climb angle: 9°
 Motor: 0.75Hp high torque motor (2 Pole)
 Batteries: 2 x 50 Amp minimum (supplied)
 Charger: 8 Amp
 Tyres (Pneumatic): Ø260mm - Front + Rear
 Product Weight: 99kg - with batteries
 Maximum User Weight: 136kg - on flat ground



0° - 136kg
 6° - 113kg
 9° - 91kg