







FEATURES

- Front and rear suspension
- Scooter pulls apart for storage and transport
- Comfortable mid back contoured seat
- Fold down adjustable steering tiller
- Straight handle with finger operated controls
- Rear wheel drive
- Electromagnetic brake
- LED headlight
- Battery power indicator
- Removable rotating / height adjustable seat
- Lift up & adjustable width arm rests
- Front bumper
- Removable shopping basket
- Interchangeable colour panels

SPECIFICATIONS

Overall: 1150mm (L) x 530mm (W) x 975mm (H)

Turning Circle: 1205mm

Speed: 6km/h

Maximum range: Up to 18km

Safe climb angle: 6° Maximum climb angle: 6°

Motor: 0.3Hp high torque motor (2 Pole) Batteries: 2 x 14 Amp minimum (supplied)

Charger: 2 Amp

Tyres (Solid black): Ø200 x 50mm / 8" - Front + Rear Maximum User Weight: 113kg - on flat ground Overall Product Weight: 46kg - with battery pack

Rear Chassis Weight: 9.5kg Front Chassis Weight: 16.8kg Battery Pack Weight: 10.5kg

Seat Weight: 8.3kg

Shopping Basket Weight: 0.9kg

shoprider.com.au

Please contact **SHOPRIDER** for more information on product