SHOPRIDER® MOBILITY SCOTERS & POWERCHAIRS

ROGW 12

BLUE

FEATURES - High back contoured se

889WBK

- High back contoured seat with headrest
- Wide seat option available
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Emergency hand brake
- Twin headlights (LED)
- Rear lights with brake-lights
- Full LCD display
- Automotive style Indicators / Hazard lights
- Large battery capacity 2 x 100Amps XT100 Premium
- Battery power indicator
- USB power outlet
- Sliding seat
- Lift up & adjustable width arm rests
- Removable shopping basket



SPECIFICATIONS

Overall: 1740mm (L) x 750mm (W) x 1310mm (H) Std. Seat: 560mm (W) x 540mm (D) - (22" x 21") Wide Seat: 820mm (W) x 450mm (D) - (32" x 18")

Turning Circle: 2705mm

Speed: 10km/h

Maximum range: Up to 50km

Safe climb angle: 9° Maximum climb angle: 9°

Motor: 2Hp high torque motor (4 Pole)

Batteries: 2 x 12V 100Ah

Charger: 8 Amp

Tyres: **Ø410mm - Front + Rear**


Product Weight: 186kg

Maximum User Weight: 350kg - on flat ground

P: **08 9248 4180**

E: orders@for-de.com.au

REDGUM Brand - 1 Business Way, Malaga WA 6090 Part of the For-De Group Pty Ltd



WIDE SEAT OPTION AVAILABLE

shoprider.com.au

FULL LCD

DISPLAY