



EZ FOLD-N-GO WALKER™



PRODUCT #4300



Folds easily for travel

COLORS



#4300-BW



#4300-RR

Weights only 8 lbs,
Supports up to 400 lbs!



#4300-CB



ACCESSORIES (SOLD SEPARATELY)



REAR GLIDES
Product # 4302



ORGANIZER POUCH
Product # 4303



SKI GLIDE SOCK COVERS
(SET OF 2) Product # 4307



REPLACEMENT WHEELS
(SET OF 2) Product # 4301

FEATURES



HANDLE HEIGHT ADJUSTMENT:
32" – 38.5"



WEIGHS ONLY 8 LBS,
SUPPORTS UP TO 400 LBS



LOCKING SWIVEL WHEELS

SPECIFICATIONS

- Weight: 8 lbs
- Collapsed Dimension: 36" H x 5.25" W x 8" D
- Wheels: 6" diameter wheels
- Case Dimensions: 37.25" L x 9.5" W x 19.25" H, 28 lbs
- Open Depth: 21.5"
- Open Width: 25.5"
- Open Handle Width: 18"
- Package Dimensions: 36.25" L x 8.5" W x 6" H, 8 lbs
- Case Dimensions: 37.25" L x 9.5" W x 19.25" H, 28 lbs
- Recommended body height range: 5'6" – 6'3"