

# Trust Care Indoor Walker

## Let's Go Indoor Walker

- Can be used in narrow spaces like the bathroom and toilet.
- Doubles as a kitchen trolley for carrying meals and coffee.
- Comes with a textile basket that can be placed under tray for storage.
- Can be operated with one hand.
- Brake may be used to park rollator to become a dinner or card tray.
- It's light, smooth and very easy to operate.
- 6" wheels, ideal for indoor use.



### Key Specifications

Handle Height	78-90cm
Tray Height	68cm
Length	63cm
Width	54cm
Net Weight	6.5kg
Max User Weight	100kg