

	Aspire Vogue Indoor Walker	Days Traymobile with Brakes	Peak Care Ellipse Indoor Walker	Trust Care Let's Go Indoor Walker																																										
Image:																																														
Description:	<p>A high-quality indoor walker with one plastic tray, a comfortable, height-adjustable handle and a one-handed brake. Ideal for users who may need to push the walker with one hand.</p> <p>Available with champagne or red frame.</p>	<p>An indoor walker with height-adjustable handles with two-way hand brakes, two trays, and four wheels. This model would suit someone who usually uses a seat walker or rollator.</p> <p>Available in beige or white.</p>	<p>A high-quality indoor walker with one plastic tray, a comfortable, height-adjustable handle and a one-handed brake. Ideal for users who may need to push the walker with one hand.</p> <p>Available with white or dark blue frame.</p>	<p>A high-quality indoor walker with one plastic tray, a comfortable, height-adjustable handle and a one-handed brake. Ideal for users who may need to push the walker with one hand.</p> <p>Available with beech or red frame.</p>																																										
Key Features:	<ul style="list-style-type: none"> • One plastic, removable tray • Foldable frame for easy storage • One foam-covered handle bar • Braking system can be used with one hand if needed • Parking brake • Comes with bag that goes under tray for extra storage • Two swivel castors at the front to easily move around corners 	<ul style="list-style-type: none"> • Two trays clip on to the frame which increases their stability • Two-way hand brakes • Handles are height-adjustable • Two swivel casters at the front to easily move around corners 	<ul style="list-style-type: none"> • One plastic, removable tray • Foldable frame for easy storage • One foam-covered handle bar • Braking system can be used with one hand if needed • Parking brake • Comes with vinyl bag that goes under tray for extra storage • Two swivel castors at the front to easily move around corners 	<ul style="list-style-type: none"> • One plastic, removable tray • Foldable frame for easy storage • One foam-covered handle bar • Braking system is easily operated and can be used with one hand if needed • Parking brake • Comes with bag that goes under tray for extra storage • Two swivel castors at the front to easily move around corners 																																										
Specifications	<table border="0"> <tr> <td><u>Handle Height</u></td> <td><u>Tray Height</u></td> <td><u>SWL</u></td> </tr> <tr> <td>86-97cm</td> <td>Not stated.</td> <td>110kg</td> </tr> <tr> <td></td> <td></td> <td><u>Tray SWL</u></td> </tr> <tr> <td></td> <td></td> <td>5kg</td> </tr> </table>	<u>Handle Height</u>	<u>Tray Height</u>	<u>SWL</u>	86-97cm	Not stated.	110kg			<u>Tray SWL</u>			5kg	<table border="0"> <tr> <td><u>Handle Height</u></td> <td><u>Tray Height</u></td> <td><u>SWL</u></td> </tr> <tr> <td>75 – 90.5cm</td> <td>Tray 1: 39cm Tray 2: 69cm</td> <td>150kg</td> </tr> </table>	<u>Handle Height</u>	<u>Tray Height</u>	<u>SWL</u>	75 – 90.5cm	Tray 1: 39cm Tray 2: 69cm	150kg	<table border="0"> <tr> <td><u>Handle Height</u></td> <td><u>Tray Height</u></td> <td><u>SWL</u></td> </tr> <tr> <td>83.5 – 94.5cm</td> <td>68cm</td> <td>100kg</td> </tr> <tr> <td></td> <td></td> <td><u>Tray SWL</u></td> </tr> <tr> <td></td> <td></td> <td>5kg</td> </tr> </table>	<u>Handle Height</u>	<u>Tray Height</u>	<u>SWL</u>	83.5 – 94.5cm	68cm	100kg			<u>Tray SWL</u>			5kg	<table border="0"> <tr> <td><u>Handle Height</u></td> <td><u>Tray Height</u></td> <td><u>SWL</u></td> </tr> <tr> <td>79 – 90cm</td> <td>68cm</td> <td>100kg</td> </tr> <tr> <td></td> <td></td> <td><u>Tray SWL</u></td> </tr> <tr> <td></td> <td></td> <td>5kg</td> </tr> </table>	<u>Handle Height</u>	<u>Tray Height</u>	<u>SWL</u>	79 – 90cm	68cm	100kg			<u>Tray SWL</u>			5kg
<u>Handle Height</u>	<u>Tray Height</u>	<u>SWL</u>																																												
86-97cm	Not stated.	110kg																																												
		<u>Tray SWL</u>																																												
		5kg																																												
<u>Handle Height</u>	<u>Tray Height</u>	<u>SWL</u>																																												
75 – 90.5cm	Tray 1: 39cm Tray 2: 69cm	150kg																																												
<u>Handle Height</u>	<u>Tray Height</u>	<u>SWL</u>																																												
83.5 – 94.5cm	68cm	100kg																																												
		<u>Tray SWL</u>																																												
		5kg																																												
<u>Handle Height</u>	<u>Tray Height</u>	<u>SWL</u>																																												
79 – 90cm	68cm	100kg																																												
		<u>Tray SWL</u>																																												
		5kg																																												