

| | Peak Care Ellipse Lite | Aspire Vogue Lightweight | Airgo Excursion | Airgo Fusion 2-in-1 Walker/Wheelchair | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|---------------|---|--|---|---|---------|-----|-------|------|-----------|-----|-------|-----|--------|-----------|-------|-------|-------|--------|------------|-------|-------|--|-------|-------------|---------------|---------|-----|-------|------|-----------|-------|-------|---|-------|-------------|---------------|---------|-----|-----|------|-----------|--------|-------|-----|------|-----------|--------|-------|-----|------|------------|--------|-------|-----|------|-----------|-------|-------|--|-------|-------------|---------------|---------|-----|--------|------|-----------|--------|-------|
| Image: | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Description: | <p>A fixed seat-height seat walker, ideal for those who travel/have limited storage space available.</p> <p>Available in red (small and medium) and black (medium and large)</p> | <p>Side-folding seat walker with a clean and crisp frame design.</p> <p>Available in red and champagne.</p> | <p>A fixed seat-height seat walker, ideal for those who travel/have limited storage space available.</p> <p>Available in at least colours from Pacific Blue, Pearl Black, Cranberry and Lime, depending upon the model.</p> | <p>This model easily converts from a seat walker to a wheelchair when required. Ideal for those who occasionally require someone to push them. Has a longer wheel base for extra stability when used as a wheelchair.</p> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Key Features: | <ul style="list-style-type: none"> Lightweight aluminium frame (approximately 8kg) Nylon bag at front Two-way hand brakes for walking and sitting. Adjustable handle height Cane holder included Frame folds sideways for easier storage and transport Highly manoeuvrable | <ul style="list-style-type: none"> Lightweight aluminium frame. Removable storage bag that holds up to 5kg. Cane holder. Two-way hand brakes for walking and sitting. Adjustable handle height. | <ul style="list-style-type: none"> Lightweight aluminium frame (approximately 7kg) Padded seat and backrest Nylon bag at front Two-way hand brakes Two-way hand brakes for walking and sitting. Frame folds sideways for easier storage and transport The XWD model is ideal for heavier users – 50cm wide and 181kg SWL | <ul style="list-style-type: none"> Lightweight aluminium frame (approximately 7.5kg without accessories) Available in blue only Padded seat and backrest Backrest flips over when converting from seat walker to wheelchair Nylon bag at front Swing away and removable foot plates Two-way hand brakes Adjustable handle height Frame folds sideways for easier storage and transport | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Models: | <table border="1"> <thead> <tr> <th>Model</th> <th>Seat Height</th> <th>Handle Height</th> <th>Net Wgt</th> <th>SWL</th> </tr> </thead> <tbody> <tr> <td>Small</td> <td>49cm</td> <td>821-945mm</td> <td>8kg</td> <td>136kg</td> </tr> <tr> <td>Med</td> <td>52.3cm</td> <td>840-978mm</td> <td>8.2kg</td> <td>136kg</td> </tr> <tr> <td>Large</td> <td>59.5cm</td> <td>931-1055mm</td> <td>8.7kg</td> <td>136kg</td> </tr> </tbody> </table> | Model | Seat Height | Handle Height | Net Wgt | SWL | Small | 49cm | 821-945mm | 8kg | 136kg | Med | 52.3cm | 840-978mm | 8.2kg | 136kg | Large | 59.5cm | 931-1055mm | 8.7kg | 136kg | <table border="1"> <thead> <tr> <th>Model</th> <th>Seat Height</th> <th>Handle Height</th> <th>Net Wgt</th> <th>SWL</th> </tr> </thead> <tbody> <tr> <td>(Med)</td> <td>54cm</td> <td>785-920mm</td> <td>6.4kg</td> <td>150kg</td> </tr> </tbody> </table> | Model | Seat Height | Handle Height | Net Wgt | SWL | (Med) | 54cm | 785-920mm | 6.4kg | 150kg | <table border="1"> <thead> <tr> <th>Model</th> <th>Seat Height</th> <th>Handle Height</th> <th>Net Wgt</th> <th>SWL</th> </tr> </thead> <tbody> <tr> <td>X18</td> <td>46cm</td> <td>750-890mm</td> <td>6.98kg</td> <td>136kg</td> </tr> <tr> <td>X20</td> <td>51cm</td> <td>770-940mm</td> <td>6.98kg</td> <td>136kg</td> </tr> <tr> <td>X23</td> <td>58cm</td> <td>850-1000mm</td> <td>7.75kg</td> <td>136kg</td> </tr> <tr> <td>XWD</td> <td>51cm</td> <td>760-930mm</td> <td>7.4kg</td> <td>181kg</td> </tr> </tbody> </table> | Model | Seat Height | Handle Height | Net Wgt | SWL | X18 | 46cm | 750-890mm | 6.98kg | 136kg | X20 | 51cm | 770-940mm | 6.98kg | 136kg | X23 | 58cm | 850-1000mm | 7.75kg | 136kg | XWD | 51cm | 760-930mm | 7.4kg | 181kg | <table border="1"> <thead> <tr> <th>Model</th> <th>Seat Height</th> <th>Handle Height</th> <th>Net Wgt</th> <th>SWL</th> </tr> </thead> <tbody> <tr> <td>Fusion</td> <td>51cm</td> <td>760-930mm</td> <td>7.76kg</td> <td>136kg</td> </tr> </tbody> </table> | Model | Seat Height | Handle Height | Net Wgt | SWL | Fusion | 51cm | 760-930mm | 7.76kg | 136kg |
| Model | Seat Height | Handle Height | Net Wgt | SWL | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Small | 49cm | 821-945mm | 8kg | 136kg | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Med | 52.3cm | 840-978mm | 8.2kg | 136kg | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Large | 59.5cm | 931-1055mm | 8.7kg | 136kg | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Model | Seat Height | Handle Height | Net Wgt | SWL | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| (Med) | 54cm | 785-920mm | 6.4kg | 150kg | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Model | Seat Height | Handle Height | Net Wgt | SWL | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| X18 | 46cm | 750-890mm | 6.98kg | 136kg | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| X20 | 51cm | 770-940mm | 6.98kg | 136kg | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| X23 | 58cm | 850-1000mm | 7.75kg | 136kg | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| XWD | 51cm | 760-930mm | 7.4kg | 181kg | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Model | Seat Height | Handle Height | Net Wgt | SWL | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Fusion | 51cm | 760-930mm | 7.76kg | 136kg | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |