





	Redgum Seat Walker Standard	Redgum Seat Walker Adjustable Seat Height	Redgum Seat Walker Push-Down Brakes	Mack Seat Walkers Heavy Duty																																													
Image:																																																	
Description:	A seat walker with adjustable handle height, pull up/lockdown handbrakes, and with a padded seat, available in a range of seat heights to suit most users.	A seat walker with adjustable seat height, adjustable handle height, pull up/lockdown handbrakes, and with a padded seat. The seat height adjustment suits most users.	A seat walker with push-down brakes, adjustable handle height, and with a padded seat. Ideal for users with severe arthritis unable to operate hand brakes. NOTE: Brakes apply automatically when user leans on handles or sits on seat.	A seat walker with increased weight capacity and wider seat, ideal for heavier individuals NOTE: Castor wheels are positioned under the frame to reduce overall width.																																													
Key Features:	<ul style="list-style-type: none"> Aluminium powder coated frame Two-way hand brakes Adjustable handle height Large vinyl under seat storage bag Comfortable folding back rest Foldable frame for easy storage and transport 6-inch or 8-inch castors depending on model 	<ul style="list-style-type: none"> Aluminium powder coated frame Two-way hand brakes Each model has a choice of 2 seat heights Adjustable handle height Large vinyl under seat storage bag Comfortable folding back rest Foldable frame for easy storage and transport 6-inch or 8-inch castors depending on model 	<ul style="list-style-type: none"> Aluminium powder coated frame Push-down brakes (brakes applied as user leans on handles or sits on seat) Adjustable handle height Large vinyl under seat storage Comfortable folding backrest Foldable frame for easy storage and transport Twin 5-inch rear castors with rubber stoppers – these provide the braking action when user leans on walker. 	<ul style="list-style-type: none"> Heavy duty powder coated aluminium frame Built-in seat pad for stability and comfort Comfortable folding back rest Two-way hand brakes Adjustable handle height Removable storage basket in front of seat Foldable frame for easy storage and transport 8-inch castors 																																													
Specifications:	<table border="1"> <thead> <tr> <th>Seat Height</th> <th>Castors Size</th> <th>SWL</th> </tr> </thead> <tbody> <tr> <td>48cm</td> <td>6-inch</td> <td>120kg</td> </tr> <tr> <td>51cm</td> <td>8-inch</td> <td>120kg</td> </tr> <tr> <td>56cm</td> <td>6-inch</td> <td>120kg</td> </tr> <tr> <td>56cm</td> <td>8-inch</td> <td>120kg</td> </tr> <tr> <td>61cm</td> <td>8-inch</td> <td>120kg</td> </tr> </tbody> </table>	Seat Height	Castors Size	SWL	48cm	6-inch	120kg	51cm	8-inch	120kg	56cm	6-inch	120kg	56cm	8-inch	120kg	61cm	8-inch	120kg	<table border="1"> <thead> <tr> <th>Seat Height</th> <th>Castors Size</th> <th>SWL</th> </tr> </thead> <tbody> <tr> <td>47cm or 56cm</td> <td>6-inch</td> <td>120kg</td> </tr> <tr> <td>52cm or 58cm</td> <td>8-inch</td> <td>120kg</td> </tr> </tbody> </table>	Seat Height	Castors Size	SWL	47cm or 56cm	6-inch	120kg	52cm or 58cm	8-inch	120kg	<table border="1"> <thead> <tr> <th>Seat Height</th> <th>Castors Size</th> <th>SWL</th> </tr> </thead> <tbody> <tr> <td>54cm</td> <td>5-inch</td> <td>120kg</td> </tr> </tbody> </table>	Seat Height	Castors Size	SWL	54cm	5-inch	120kg	<table border="1"> <thead> <tr> <th>Model</th> <th>Seat Height</th> <th>SWL</th> </tr> </thead> <tbody> <tr> <td>Mini Mack</td> <td>50cm</td> <td>225kg</td> </tr> <tr> <td>Maxi Mack</td> <td>58cm</td> <td>225kg</td> </tr> <tr> <td>Mighty Mack (wider seat)</td> <td>55cm</td> <td>225kg</td> </tr> </tbody> </table>	Model	Seat Height	SWL	Mini Mack	50cm	225kg	Maxi Mack	58cm	225kg	Mighty Mack (wider seat)	55cm	225kg
Seat Height	Castors Size	SWL																																															
48cm	6-inch	120kg																																															
51cm	8-inch	120kg																																															
56cm	6-inch	120kg																																															
56cm	8-inch	120kg																																															
61cm	8-inch	120kg																																															
Seat Height	Castors Size	SWL																																															
47cm or 56cm	6-inch	120kg																																															
52cm or 58cm	8-inch	120kg																																															
Seat Height	Castors Size	SWL																																															
54cm	5-inch	120kg																																															
Model	Seat Height	SWL																																															
Mini Mack	50cm	225kg																																															
Maxi Mack	58cm	225kg																																															
Mighty Mack (wider seat)	55cm	225kg																																															